



Pacer

Gait training has come a long way since we introduced our first model in 1977, and the Pacer continues to evolve and improve to meet our customers' needs. Over the years your suggestions have helped us keep our design at the cutting edge.

The Pacer, because it can be configured in so many ways, can accommodate clients of every ability, in any environment. It's the versatility and dependability you've come to expect from Rifton.



Find letters of medical necessity and informative articles at: www.rifton.com/pacer

students have a wide range of motor issues. Many need a gait trainer that allows them more freedom of movement. Sometimes with our old equipment, it seems like they're fighting against the device, but Rifton's dynamic Pacer allows and even shapes the natural movement that occurs with walking. It's fabulous.

Nikki Cornell, MPT Clovis Unified School District, CA

The Pacer (dynamic shown here) is easily configured to meet your client's every need. Choose between dynamic and standard upper frames and pair them with standard bases.



Pacer features


The six members of the Pacer family











Perfect positioning made easy

The single-column upper frame

features a locking gas spring that adjusts to any height on the indexed column.





Combining the fluid height adjustments of the MPS and upper frame brings the seat low enough for a transfer from the seated position.



Once the transfer is completed, caregivers raise the upper frame to bring the client into an upright position.



The multi-position saddle (MPS)

is the ultimate pelvic positioning aid. Adjustable in five directions to accommodate each client's posture, the MPS easily attaches to both the dynamic and standard upper frames.





Hip corral depth and height adjustments move the hip corral to provide cushioning and support for the pelvis at the level of the ASIS.

Adjustment range: Large: depth 3", height 3" Med/Small: depth 3", height 2"



Saddle depth adjusts to optimally position the pelvis in relation to the upper body to encourage weight-bearing and forward motion during ambulation.

Adjustment range: Large: 5", Med/Small: 4"



Saddle height adjusts to fine-tune the amount of weight-bearing support during gait.

Adjustment range:

Large: 81/2", Med: 81/2", Small: 51/2"

Saddle angle is adjustable to position the pelvis in anterior or posterior tilt for the most effective gait pattern.

Adjustment range:

Large: 15° forward and 15° back Med/Small: 7° forward and 7° back





Using the height adjustment on the MPS, the caregivers raise the client's pelvis to the optimal height for weight-bearing and step-taking.



With the tilt adjustment on the MPS, the caregivers position the client's pelvis for the most effective gait pattern.

Infinite control for effective movement

The Pacer casters' five functions let you fine-tune the

Pacer's movement as your clients gain control.

Standard base casters



Utility base front caster

Utility base rear wheel





Utility base

Easy rolling over thresholds or thick carpet (or even outdoors on gravel or grass) is possible with the large-wheeled utility base, letting you go places you never could before.





Transportable

The detachable upper frame makes it easy to stow and transport the Pacer. Separable components allow you to bring the right base for each client and every setting.



Go minimal

The Pacer chest pad and handlebars form a lightweight configuration perfect for higher functioning users. Available for small, medium and large Pacer sizes, these components allow easy transfers and quick adjustment.







Posterior position



Posterior positioning is a natural progression towards independent movement, and it gives the users improved access to their environment.

Note: The standard base is recommended for posterior positioning. When the utility base is reversed, steering is more difficult and the directional lock feature will not work.

Reverse the MPS

The MPS has been designed to allow reverse positioning of the saddle, which enables gait practice with the frame in the posterior position. For instructions on how to place the MPS in the posterior position, visit www.rifton.com/MPSposterior.



Components



Dynamic upper frame provides dynamic weight-shifting and dynamic body-weight support, encouraging natural gait patterns.



Standard upper frame is a light-weight choice for clients who do not need dynamic movement.



Standard base

Designed to be used on smooth surfaces.

Caster size: 5½"

Utility base

For indoors or out, this base is the best choice for navigating thresholds, lawns, gravel paths or chipped playgrounds. No downhill though! Always be safe.

Caster/wheel size: front casters 8", rear wheels 11½"





Arm prompts

Versatile arm prompts allow height, rotation, angle, lateral proximity and forward/backward adjustments to accommodate numerous positioning requirements.

Padded surface area Adjustment range small 4½"x 8" (elbow to fist): small 9'

small 4½"x 8" (elbow to fist): small 9"-12½' large 5"x 9" large 10½"-15"





Arm platforms

Arm platforms offer a simple alternative for clients who need less arm positioning. They attach to the adjustable arm prompt base and are made of closed-cell molded foam for padded comfort and easy cleaning.

Padded surface dimensions: $13\frac{1}{2}$ "x $5\frac{1}{2}$ "





Hand brakes

The hand brakes combine both running and parking brake functions in a simple leveraction control. The brakes attach to the arm prompts and can be quickly removed when not needed.

Note: Since retrofitting is complex, we recommend ordering brakes with the original purchase.





Hand loops

For less involved clients, hand loops are a simpler option for upper extremity support.





Chest prompt

The padded chest prompt holds the client's torso securely at the angle of your choice. It is available in three sizes.

Range of circumference	Pad heigh
small 14"-28"	small 4½"
med 22"-40"	med 61/2"
large 28"-50"	large 8"



Components (continued)



Chest pad

The chest pad is ideal for less involved users. With its single support strap and double angle adjustment it allows easy transfers and quick positioning.

Range of circumference	Pad height
small 14"-28"	small 4"
med 22"-40"	med 5"
large 28"-50"	large 6"





Handlebars

For users who need less arm support, the handlebars attach quickly and adjust easily. They pair perfectly with the simple chest pad.

Height from pivot to top: small: 5", large: 7"





Multi-position saddle

The multi-position saddle (MPS) enables precise positioning of the pelvis. Adjustable in five directions to accommodate each client's posture, the MPS attaches to both the dynamic and standard upper frames. A seat cover is also available to decrease surface friction.

Dimensions: see dimensions chart





Pelvic support

For clients who need less abduction, the pelvic support is a softer alternative to the hip positioner. It provides weight-bearing assistance and freedom of movement while clients build strength and stamina.

Dimensions: see dimensions chart





Hip positioner

The hip positioner enables natural pelvic movement and lower extremity weight-shifting, and provides good abduction. A pad is available for added comfort.

Dimensions: see dimensions chart







Thigh prompts


Thigh prompts swing with the user's stride and prevent leg scissoring. They are adjustable in the medial and lateral directions and prevent the user from turning within the Pacer frame.

Length (from clamp to strap) small 5"-8" large 5"-11"

Max. leg circumference

small 15" large 20"





Ankle prompts

Ankle prompts have comfortably padded ankle straps and spring adjusters to control stride length and placement. Adjust the strap length to determine foot separation and prevent scissoring.





Handholds

The ergonomic handholds are heightadjustable, and can be placed anywhere along the top bar of the frame. Handholds can be used by caregiver or client.





Communication tray

The communication tray is great for holding tablets and other mobile devices.

Inside dimensions: 8½"x 11" Max. working load: 10 lb





Attendant guide bar

The guide bar is made of lightweight, high strength aluminum alloy. It attaches to the front tube of the standard or utility base with a quick-release clamp, and enables the caregiver to guide the Pacer from in front or behind.





Accessories tote

A tough and attractive tote bag big enough (14"x 16"x 9") to keep all your Pacer components together when not in use. Made of black ballistic nylon, here is a bag that will hold up to whatever you throw in it and still look great.



Pacer dimensions


Key user d	imensions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL			
Elbow height			151/2-201/2	181/2-271/2	24-35	32-47	34–49			
	Key user dimension: Measure the vertical d while the user is stand	istance from the bent elbov	v to the floor							
八	Choose the size that allows for growth.									
Item dimens	ions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL			
Floor to top of	arm prompt		151/2-201/2	181/2-271/2	24-35	32-47	34-49			
Floor to top of	MPS (not tilted)*		N/A	10 – 20	91/2 - 231/2	17 – 35	191/2-371/2			
Overall width		Standard base:	201/2	23	26	28	311/2			
		Utility base:	N/A	N/A	30	32	N/A			
Overall length		Standard base:	221/2	26	30	36	39¾			
		Utility base:	N/A	N/A	36	401/2	N/A			
Frame height		Standard & utility base:	121/2	16-21	21½-27	281/2-381/2	31-41			
Base height (w	thout upper frame)	Standard & utility base::	N/A	11	141/2	16	181/2			
Frame weight (b)	Standard base:	101/4	11½	15	151/2	221/2			
		Utility base:	N/A	N/A	181/2	191/4	N/A			
		Standard upper:	N/A	6	71/2	91/2	11			
		Dynamic upper:	N/A	71/2	11½	131/2	15			
Dynamic upper	movement	Vertical:	N/A	11/4	2	2	2			
		Horizontal:	N/A	11/2	2	2	2			
Maximum work	ring load (lb)		50	75	150	200	250			

^{*}This measurement is taken at the lowest frame and MPS height, and the highest frame and MPS height.

20

12

Pelvic support



16


9

12

Hip positioner

Dimensions (inches)—without pad				
Rear	. <u>.</u> <u>.</u> <u>.</u> <u>.</u> <u>.</u> <u>.</u> <u>.</u> <u>.</u> <u>.</u> .	Front		
	K513 small	K533 large		
Crotch length	21/2	31/2		
6 Width	11/4	13/4		
G Length	9	11		

MPS saddle





All five sizes are available in five colors: red, blue, gray, lime and pink.



G

2024 order form

To order or get a quote:
Go online: https://www.rifton.com/pacer-sizing
Email: sales@rifton.com
Call: 800.571.8198

Use dimension chart to select appropriate size.

				K610 mini		K620 small		K630 medium	K640 large	K650 X	L
* Color		Red		Red		Red 🗆	Red [Red			
				Blue		Blue		Blue	Blue	Blue	
			CHOOSE	Gray		Gray		Gray 🔲	Gray [] Gray	
			SE	Lime		Lime		Lime	Lime] Lime	
				Pink		Pink		Pink 🔲	Pink [] Pink	
* Up	per frame	웃	Dynamic	N/A		K621		K631	K641] K651	
		CHOOSE	Standard	N/A		K622		K632	K642] K652	
	Standard base	운	With odometer	K613		K623		K633	K643] K653	
*		CH00SE	Without odometer	K614		K624		K634 🔲	K644] K654	
Base	Utility base	운	With odometer	N/A		N/A		K635	K645] N/A	
	(larger wheels)	CHOOSE	Without odometer	N/A		N/A		K636	K646 [] N/A	
Arm	supports (pair)		Small arm prompts	K512		K512		K512	K512] K512	
		오	Large arm prompts	N/A		N/A		K532	K532] K532	
		CHOOSE	Arm platforms	N/A		N/A		K697	K697	-	
		•	Arm platforms w/ handgrips	N/A		N/A		K696 🗆	K696		
Hand	Ibrake for utility base (requires a	rm sup		N/A		N/A		K687	K688		
	llebars (pair)	•	Small	K625		K625		K625	K625		
	4. /	CHOOSE	Large	N/A		N/A		K675	K675	-	
Hand	l loops (pair)		1 3	K579		K579		K579	K579		
	t prompt		Small	K514		K514		K514	N/A	N/A	
Chest pad		요		N/A		K524		K524	K524		
		3SC		N/A		N/A		K534	K534		
		N/A		K629		K639	K649				
	i-position saddle		Small	N/A		K672		N/A	N/A	N/A	
		CHOOSE	Medium	N/A		N/A		K673	N/A	N/A	
		3SC	Large	N/A		N/A		K674	K674		
MPS seat cover		Small/Medium	N/A		K627		K627	N/A	N/A		
		Large	N/A		N/A		K628	K628			
Pelvic support (includes handholds)			Small	K556		K556		K556	N/A	N/A	
		CHOOSE	Medium	N/A		K557		K557	K557		
		SE	Large	N/A		N/A		K558	K558		
Hip r	positioner (includes handholds)	요	Small	K513		K513	П	K513	K513		
inp positioner (includes handholds)		IOOSE	Large	N/A		N/A		K533	K533		
Hip positioner pad			Small	K541		K541		K541 🗆	K541		
		CHOOSE	Large	N/A		N/A		K542	K542		
Thigh prompts (pair)			Small	N/A		K515		K515	K515		
9	·· b· ambra /ban/	CHOOSE	Large	N/A		N/A		K535	K535		
Ankl	e prompts for standard base (pa	- "	25.90	K681		K682		K683			
Ankle prompts for utility base (pair)		N/A		N/A		K693					
Additional handholds (pair) Communication tray			K519		K519		K519				
			K513		K698		K698				
	ndant guide bar			K589		K679		K689			
Rifton accessories tote			1,303	ш	1073	П	1,000	1,000	7 13//		

^{*} Required

Adult gait training options

XL Pacer

The **XL Pacer** has dynamic body support, locking gas-spring height adjustment and a multi-position saddle.

E-Pacer

The **E-Pacer** combines an electric lift column with gait training features from the Pacer.

TRAM

The **TRAM** has an electric lift column and a compact frame that is easy to maneuver in tight spaces.





Product comparison chart







XL Pacer

E-Pacer

TRAM

Weight limit	250 lb	350 lb	350 lb
Width of base (to determine min. door width)	31½"	311/2"	27½"
Adjustment range (floor to top of arm prompt)	34"-49"	24"-61"	23"-60"
Zero-lift transfers (powered height adjustment)		√	√
Bluetooth-equipped scale		✓	√
Expandable base frame		√	√
Low-base option (to fit under tight spaces)			√
Collapsible frame	√		
Dynamic weight shifting	√		
Thigh prompts	√		
Multi-position saddle (MPS)	√		
Odometer	√	√	
Ankle prompts	√	✓	
Four-function casters	√	√	
Hip positioner	√	√	√
Pelvic support	✓	✓	√
Arm prompts	✓	√	√