



Stander Prone/Supine/Multi-Position

It's finally here – the Stander that has it all. We're thrilled to introduce the first two sizes of a total redesign of our Supine and Prone Stander families. This new Stander can be configured as supine, prone or multi-position.

All models allow hip abduction up to 30° to promote hip joint development. Independent adjustments for hip and knee flexion and leg length accommodate contractures and discrepancies. A broad range of positioning components allow users to remain comfortably upright and active for longer. The frame can be raised and lowered to position a child at eye level with peers, or in easy reach of activities.

All this function is contained within an elegant and approachable design that centers the child while opening worlds of possibility.



Find important details and studies on the therapeutic benefits of standing at: www.rifton.com/standers

ability to fully customize joint positioning through the lower extremities for maximal comfort and alignment is a game changer! I have a student who only tolerated five to ten minutes of static standing in a less adjustable stander. When we tried the new Rifton Stander the student was happy and comfortable for nearly an hour, and we only got him out because we ran out of time.

Linda K. Miller, PT, DPT **MOVE International Trainer** Prospect Center, NY

•••••• Size 1 Supine Stander configuration (S410)

Stander features



Size 2 Supine Stander



Multi-Position Stander



The Stander can be purchased in the prone, supine or multi-position configuration. If the multi-position configuration is ordered, the conversion between prone and supine is quick, easy and tool free.

Quick and easy conversion between prone and supine



- Add/Remove the upper trunk support.
- **2 Rotate** the tray.
- Swap the supine hip strap and prone pelvic harness.
- 4 Add/Remove the supine knee straps.
- **5** Turn the sandals 180°.



Supine

Quick, precise positioning



Advanced legs have independent leg length, hip angle and knee angle adjustments.



Up to 30° of abduction.



Tilt adjustments from 0°-90°.



Hip abduction

Independent leg abduction, up to 30° in each leg, is now standard. It's what clinicians everywhere asked us for to promote healthy hip development.





The Size 2 Stander advanced legs have independent hip and knee angle adjustment to accommodate contractures. Separate upper and lower leg length adjustments ensure that users with leg length discrepancies are comfortably supported.







Move between eye-level and activity level within seconds.

- ▲ When the frame is lowered, users can interact at eye-level with their peers.
- With the frame raised, users can access activities and environments.




Quick, safe, tool-free transfers



Supine







Components



Knee cuffs (Prone only)

The knee cuffs provide lateral stabilization and rotational control and prevent knee hyperextension.





Head laterals (Supine only)

The head laterals mount on the headrest to provide lateral head support. They adjust independently in height and width.





Butterfly harness (Supine only)

The butterfly harness can be used instead of the trunk strap for greater upper trunk stabilization.





Hip strap (Supine only)

This strap supports the user at the level of the hips. The rotation control version provides rotational stabilization.





Split knee supports

The soft knee straps form around the knee and have a cutout for the patella. For more rigid support, choose the adjustable split knee supports which swing away for easy loading.







Trunk strap

The trunk strap supports the user's chest. The rotation control version wraps around the user and provides rotational stabilization.





Leg types

Standard legs allow for hip abduction, and independent leg length adjustments. Advanced legs are optional on Size 2 and include hip and knee flexion in both prone and supine configurations.





Tray

The tray can be used in prone and supine configuration. Its position and angle adjust easily, and it swings away for transfer.





Additional trunk laterals

The laterals attach at chest or hip level and swing away for transfer. The pads rotate for individual height adjustment.





Medial/lateral thigh supports

The medial/lateral thigh supports provide abduction or adduction. The mounting bracket is off-center, so the pad has two possible height positions depending which side it is installed on.





Sandal lift kit

The smallest users may require the \$409 Sandal lift kit for best positioning.



Stander dimensions

User dimensions (inches)	S410 Size 1	S420 Size 2						
Height	28-42*	38-53						
Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth. *The smallest users may require the S409 Sandal lift kit for best positioning. Important: User's weight must not exceed the maximum working load.								
Item dimensions (inches)	S410 Size 1	S420 Size 2						
Base length x width	32½ x 21½	39 x 27						
Height of horizontal body support above floor	22-28	221/2-361/2						
Sandal to knee height	63/4-11	10-151/4						
Knee to hip height	61/4-101/2	81/2-131/2						
Sandal to top of headrest	31½-43	41–55						
Width between trunk laterals	51/2-10	7-111/2						
Width between medial/lateral thigh supports	21/2-4	31/4-5						
Basic item weight PRONE (lb)	47	72						
Basic item weight SUPINE (lb)	52	78						
Max. working load (lb)	75	100						







Size 2 range of adjustment



2024 order form

Use dimension chart to select appropriate size.



					S410 size 1		S420 size 2	
* Frame					S411		S421	
* Color Google			Blue		Blue			
				ÖSE	Red		Red	
ж СНО	Prone (with laterals, pelvic harness)			S412		S422		
	HCPCS code: E0638	Prone knee cuffs (pair)			S485		S485	
	Supine (with upper trunk				S413		S423	
	support, laterals)	Head laterals (pair)			S497		S497	
	HCPCS code: E0638	Butterfly harness			S456		S457	
OSE		* Hip strap	CHOOSE	Standard	S472		S473	
CHOOSE CONFIGURATION			OSE	Rotation control	S476		S477	
NFIC		* Knee Supports	CHOOSE	Soft straps (pair)	S401		S402	
jUR/			OSE	Split supports (pair)	N/A		S487	
OITA	Multi-Position (with upper				S414		S424	
Ž	trunk support, laterals, pelvic harness)	Head laterals	lead laterals		S497		S497	
	HCPCS code: E0641	Butterfly harness	Butterfly harness		S456		S457	
	* Hip s	* Hip strap	CHOOSE	Standard	S472		S473	
			SE	Rotation control	S476		S477	
		* Knee Supports	СНООЅЕ	Soft straps (pair)	S401		S402	
			OSE	Split supports (pair)	N/A		S487	
		Prone knee cuffs (pair)		S485		S485		
* Trunk strap Standard Rotation control				S462		S463		
			OSE	Rotation control	S466		S467	
* Leg type			CHOOSE	Standard	S415		S425	
			3SO	Advanced	N/A		S426	
Tray				S491		S491		
Additional trunk laterals (pair)				S494		S494		
Medial/lateral thigh supports (pair)				S482		S482		
Additional medial/lateral thigh supports (pair)				S482		S482		
S1 Sandal lift kit				S409		NA		

The HCPCS codes indicated above are suggestions only, based on knowledge of our products and the HCPCS definitions established by CMS. It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. Actual coding and coverage varies by payer, some may accept or require alternative HCPCS codes, including miscellaneous codes, to ensure access for their beneficiaries.