

Positioning Checklist for the Rifton Pacer

Individual's name _____

Staff name _____

Date _____

Pacer Size

- mini
- small
- med
- large
- XL

Frame Orientation

- forward facing/Standard
- backward facing/Reverse

Frame Height

Distance between white adjustment lever and top of black molding on Pacer upright:

Chest Prompt



Size small medium large

Clamp Position (post holder)

right side of Pacer frame inside frame outside frame

left side of Pacer frame inside frame outside frame

Height Adjustment

Specify number of notches visible above clamp _____

Width Adjustment

Distance between upright post and chest prompt _____

Fit of chest prompt

snug to trunk loose to trunk

Is the chest prompt flipped? (Use the imprinted size indication on the black plastic prompt holder as a guide. If it is upside down, the chest prompt is flipped.)

yes no

Orientation of chest prompt

vertical tilted forward tilted backward

When using SoloLift to place individual in Pacer:

vest removed while in Pacer vest is on while in Pacer

Left Arm Prompt/Hand Loop



hand loops



arm prompts



NOTE: Describe all positioning as if you are standing behind the individual

Size small large

Hand loop being used? (in place of arm prompt)

yes no

Orientation of hand loop

vertical tilted forward tilted backward

Clamp placement (L-bar holder)

front bar side bar

Clamp position (L-bar holder)

inside frame outside frame

Distance of clamp from end of top bar: _____

L-bar post position

points forward points backward points right points left

L-bar height

Specify number of notches visible above clamp: _____

Orientation of arm prompt

horizontal angled up angled down

Rotation position (specify as o'clock when standing behind individual in the Pacer)

Individual needs velcro strap behind elbow? (This is to prevent the arm from slipping off the prompt. If "No", assume the Velcro strap is used in standard configuration across the forearm.)

yes no

Right Arm Prompt/Hand Loop



hand loops



arm prompts



NOTE: Describe all positioning as if you are standing behind the individual

Size small large

Hand loop being used? (in place of arm prompt)

yes no

Orientation of hand loop

vertical tilted forward tilted backward

Clamp placement (L-bar holder)

front bar side bar

Clamp position (L-bar holder)

inside frame outside frame

Distance of clamp from end of top bar: _____

L-bar post position

points forward points backward points right points left

L-bar height

Specify number of notches visible above clamp: _____

Orientation of arm prompt

horizontal angled up angled down

Rotation position (specify as o'clock when standing behind individual in the Pacer)

Individual needs velcro strap behind elbow? (This is to prevent the arm from slipping off the prompt. If "No", assume the Velcro strap is used in standard configuration across the forearm.)

yes no

Pelvic Support



pelvic support



Size small medium large

Front strap attachment

front bar side bar chest prompt post

Length of front straps: _____

Clamp position (for back strap attachment)

inside frame outside frame

Distance from end of top bar: _____

Handhold height (for back strap attachment)

Specify number of notches visible above clamp: _____

Handhold position (for back strap attachment)

points forward points backward

Length of back straps: _____

..... **OR**

Hip Positioner



hip positioner

Size small large

Use with hip positioner pad? yes no

Front strap attachment

front bar side bar chest prompt post

Length of front straps: _____

Clamp position (for back strap attachment)

inside frame outside frame

Distance from end of top bar: _____

Handhold height (for back strap attachment)

Specify number of notches visible above clamp: _____

Handhold position (for back strap attachment)

points forward points backward

Length of back straps: _____

Thigh Prompts



Size small large

Use on:

right thigh left thigh both

Clamp position

Right leg/side of Pacer frame

inside frame outside frame

Left leg/side of Pacer frame

inside frame outside frame

Ankle Prompts



Use on:

right ankle left ankle both

Ankle prompt placement:

Right side of Pacer frame

Distance of front spring from front caster attachment _____

Distance of back spring from back caster attachment _____

Length of strap from frame _____

Left side of Pacer frame

Distance of front spring from front caster attachment _____

Distance of back spring from back caster attachment _____

Length of strap from frame _____

Casters



Front caster settings:

no setting needed

swivel lock

set ratchet (Push down white ratchet lever in the center. This keeps wheels from rolling backwards.)

set drag (Move white dial towards the turtle imprinted on the plastic. Dial setting: _____)

Rear caster settings: (see explanations above)

no setting needed swivel lock set ratchet set drag

(Dial setting: _____)