

Positioning Checklist for the Rifton Activity Chair

Use this Positioning Checklist as a convenient way to assure optimal use of the Rifton Activity Chair for each individual. Write notes to customize your instructions for each individual.

INDIVIDUAL'S NAME _____

STAFF NAME _____

DATE _____

THIS POSITIONING CHECKLIST IS INTENDED FOR CONSISTENT USE OF ADAPTIVE EQUIPMENT FOLLOWING ASSESSMENT BY A HEALTH PROFESSIONAL. THE CHECKLIST DOES NOT REPLACE AN INITIAL CLINICAL EVALUATION FOR CORRECT SIZING, COMPONENT SELECTION, AND APPROPRIATE THERAPEUTIC POSITIONING FOR THE INDIVIDUAL. SOME COMPONENT SIZING IS DEPENDENT ON THE BASE MODEL SIZE.



Hi/Lo Base

Size

- Small
- Medium
- Large



Standard Base

Size

- Small
- Medium
- Large

Height of legs _____

Backrest Spring Action



- Locked
- Unlocked

Comments _____

(Some individuals may need to have it locked during specific activities.)

Base Spring Action (available on standard base only)



- Locked
- Unlocked

Comments _____

(Some individuals may need to have it locked during specific activities.)

Tilt in Space



Measure the angle (using the protractor on side of the seat)

Tilt forward _____

Tilt backward _____

Comments _____

(Some individuals may need to have the seat angled forward during specific activities.)

Seat Depth



Inches (cm) _____

Seat Width



Hip Guides

Inches (cm) _____

Seat Back



Height inches (cm) _____

Measure the angle (use backrest protractor)

Forward _____

Backward _____

Footrest



Footboard

Knee angle: Length of exposed metal rod on footboard gas spring _____

Footboard height: Specify number of notches visible _____

Footplate angle: Acute Neutral Obtuse



Sandals

Left foot Right foot



Wedges

Left foot Right foot

Comments _____



Ankle straps

Left foot Right foot

Armrests



Angle of armrests:

Flat Up Down

Comments _____

Tray



With handhold

Angle tray:

Flat Up Down

Arm Prompts



Left side

Clamp position

Inside Outside

L-bar support position

Points forward Points backward
 Points toward right Points toward left

L-bar height

Specify number of notches visible above clamp _____

Arm Prompt

Angle up Angle down Straight

Rotation position for left side

Specify as o'clock when sitting in chair _____

Right side

Clamp position

Inside Outside

L-bar support position

Points forward Points backward
 Points toward right Points toward left

L-bar height

Specify number of notches visible above clamp _____

Arm prompt

Angle up Angle down Straight

Rotation position for right side

Specify as o'clock when sitting in chair _____

Seat Belt



Connect in which slot

For small or medium Activity Chair

- Forward slot
- Rear slot

For large Activity Chair

- Forward slot
- Middle slot
- Rear slot

Pelvic Harness



Connect in which slot

- Forward slot
- Middle slot (large chair only)
- Rear slot

Thigh Strap



Leg Prompts



Right leg

- In Middle Out

Left leg

- In Middle Out

Abductor



- In
- Middle
- Out

Adductors



Right leg

- In Middle Out

Left leg

- In Middle Out

Headrest



Winged Contoured Flat

Height of headrest – inches (cm) _____
(Measure from seat to top of headrest)



Laterals



Width apart in inches (cm) _____

Height of right lateral – inches (cm) _____
(Measure from seat to top of support)

Height of left lateral – inches (cm) _____
(Measure from seat to top of support)

With chest strap: Yes No

Size of laterals: Large Small

Wide Chest Strap



Height – inches (cm) _____
(Measure from seat to top of chest strap)

Butterfly Harness

Standard Slim cut



Connect in which slot

- Forward slot
- Middle slot (large chair only)
- Rear slot