

2024 order form

Use dimension chart to select appropriate size.



To order or get a quote:
 Go online: <https://www.rifton.com/wave-sizing>
 Email: sales@rifton.com
 Call: 800.571.8198

Conversion kits

The tool-free conversion kit converts your small chair to medium as your child grows.

Kit without calf rest

Blue Z268 \$390
 Pink Z269 \$390

(Kit includes seat, back, cover and chest strap)



Kit with calf rest

Blue Z266 \$485
 Pink Z267 \$485


(Kit includes seat, back, cover and chest strap)



		Z210 small	Z220 medium	Z230 large	
* Chair		Z211 \$320 <input type="checkbox"/>	Z221 \$330 <input type="checkbox"/>	Z231 \$410 <input type="checkbox"/>	
* Fabric cover	CHOOSE	Blue	Z244 \$235 <input type="checkbox"/>	Z246 \$245 <input type="checkbox"/>	Z248 \$260 <input type="checkbox"/>
		Pink	Z245 \$235 <input type="checkbox"/>	Z247 \$245 <input type="checkbox"/>	Z249 \$260 <input type="checkbox"/>
* Chest strap	CHOOSE	With lateral positioning	Z212 \$115 <input type="checkbox"/>	Z222 \$125 <input type="checkbox"/>	Z232 \$125 <input type="checkbox"/>
		Without lateral positioning	Z213 \$105 <input type="checkbox"/>	Z223 \$115 <input type="checkbox"/>	Z233 \$115 <input type="checkbox"/>
Head blocks		Z214 \$110 <input type="checkbox"/>	Z214 \$110 <input type="checkbox"/>	Z234 \$110 <input type="checkbox"/>	
Leg straps		Z215 \$105 <input type="checkbox"/>	Z225 \$115 <input type="checkbox"/>	Z235 \$125 <input type="checkbox"/>	
Calf rest (Calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216 \$145 <input type="checkbox"/>	Z226 \$160 <input type="checkbox"/>	Z236 \$185 <input type="checkbox"/>	
Tub stand		Z217 \$330 <input type="checkbox"/>	Z217 \$330 <input type="checkbox"/>	Z217 \$330 <input type="checkbox"/>	
Shower stand		Z218 \$875 <input type="checkbox"/>	Z218 \$875 <input type="checkbox"/>	Z218 \$875 <input type="checkbox"/>	
Tub transfer base (legs required, calf rest recommended)	CHOOSE	Counterclockwise into tub	Z227 \$1,035 <input type="checkbox"/>	Z227 \$1,035 <input type="checkbox"/>	Z227 \$1,035 <input type="checkbox"/>
		Clockwise into tub	Z228 \$1,035 <input type="checkbox"/>	Z228 \$1,035 <input type="checkbox"/>	Z228 \$1,035 <input type="checkbox"/>
Legs for tub transfer base	CHOOSE	Standard	Z224 \$120 <input type="checkbox"/>	Z224 \$120 <input type="checkbox"/>	Z224 \$120 <input type="checkbox"/>
		Long (tub edge over 20" high)	Z238 \$120 <input type="checkbox"/>	Z238 \$120 <input type="checkbox"/>	Z238 \$120 <input type="checkbox"/>
Additional head blocks		Z214 \$110 <input type="checkbox"/>	Z214 \$110 <input type="checkbox"/>	Z234 \$110 <input type="checkbox"/>	
Additional fabric cover		Z219 \$235 <input type="checkbox"/>	Z229 \$245 <input type="checkbox"/>	Z239 \$260 <input type="checkbox"/>	
Additional chest strap	CHOOSE	With lateral positioning	Z212 \$115 <input type="checkbox"/>	Z222 \$125 <input type="checkbox"/>	Z232 \$125 <input type="checkbox"/>
		Without lateral positioning	Z213 \$105 <input type="checkbox"/>	Z223 \$115 <input type="checkbox"/>	Z233 \$115 <input type="checkbox"/>

* Required

Wave dimensions

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40–56	50–74
 <p>Key user dimension: height Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth.</p> <p>Important: Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.</p>			
	Item dimensions (inches)	Z210 small	Z220 medium
Seat depth	10½	13	17
Seat width	15	15	17½
Back height	23	26½	32
Seat height	2–7	2–7	2–7
Seat height (incl. tub stand)	12½–18½	12½–18½	12½–18½
Seat height (incl. shower stand)	27½	27½	27½
Seat height (incl. tub transfer base) with long leg inserts	21–26 24½–29½	21–26 24½–29½	21–26 24½–29½
Length of mesh area when flat (without calf rest)	35	41	50½
Length of mesh area when flat (with calf rest)	42½	49	61
Base width	18	18	21
Base length (legs down) with legs flat	24 29½	24 29½	24 29½
Tub stand height (legs flat)	8	8	8
Tub stand length (legs flat)	39	39	39
Tub stand width (legs flat)	17½	17½	17½
Shower stand height	22	22	22
Shower stand width	27	27	27
Shower stand length	34	34	34
Tub transfer base footprint	28 x 30	28 x 30	28 x 30
Max. tub height for transfer base with long leg inserts	20 24	20 24	20 24
Maximum working load (lb)	75	150	250
Back and seat angle adjustments (degrees)		0°, 15°, 30°, 45°, 60°, 75°	



All three sizes are available in blue and pink.