

2023 order form for large

Use dimension chart to select appropriate size.



		Standard R860 large		Hi/lo R870 large		
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861	\$730 <input type="checkbox"/>	R871	\$1,390 <input type="checkbox"/>
		Adjustable w/ back spring	R862	\$845 <input type="checkbox"/>	R872	\$1,510 <input type="checkbox"/>
* Base	CHOOSE	Without spring	R873	\$445 <input type="checkbox"/>	R883	\$3,015 <input type="checkbox"/>
		With spring	R874	\$535 <input type="checkbox"/>	N/A	
* Arm supports (pair)	CHOOSE	Armrests	R868	\$260 <input type="checkbox"/>	R868	\$260 <input type="checkbox"/>
		Forearm prompts	R879	\$675 <input type="checkbox"/>	R879	\$675 <input type="checkbox"/>
* Adjustable legs	CHOOSE	Short	R808	\$70 <input type="checkbox"/>	N/A	
		Long	R894	\$85 <input type="checkbox"/>	N/A	
		Short with casters	R895	\$315 <input type="checkbox"/>	N/A	
		Long with casters	R896	\$340 <input type="checkbox"/>	N/A	
* Pads	CHOOSE	Tan	R875	\$270 <input type="checkbox"/>	R875	\$270 <input type="checkbox"/>
		Blue	R876	\$270 <input type="checkbox"/>	R876	\$270 <input type="checkbox"/>
		Red	R877	\$270 <input type="checkbox"/>	R877	\$270 <input type="checkbox"/>
		Pink	R977	\$270 <input type="checkbox"/>	R977	\$270 <input type="checkbox"/>
		Purple	R978	\$270 <input type="checkbox"/>	R978	\$270 <input type="checkbox"/>
		Green	R979	\$270 <input type="checkbox"/>	R979	\$270 <input type="checkbox"/>
Headrest	CHOOSE	Contoured	R897	\$250 <input type="checkbox"/>	R897	\$250 <input type="checkbox"/>
		Flat	R898	\$250 <input type="checkbox"/>	R898	\$250 <input type="checkbox"/>
		Adjustable winged	R881	\$275 <input type="checkbox"/>	R881	\$275 <input type="checkbox"/>
Laterals (pair)	CHOOSE	Small	R828	\$210 <input type="checkbox"/>	R828	\$210 <input type="checkbox"/>
		Large	R838	\$225 <input type="checkbox"/>	R838	\$225 <input type="checkbox"/>
		Large with chest strap	R839	\$270 <input type="checkbox"/>	R839	\$270 <input type="checkbox"/>
Push handles (recommended for all chairs with wheels)			R878	\$160 <input type="checkbox"/>	Included	
Tray (requires armrests)			R669	\$305 <input type="checkbox"/>	R669	\$305 <input type="checkbox"/>
Handhold			R649	\$30 <input type="checkbox"/>	R649	\$30 <input type="checkbox"/>
Chest strap (wide)			R807	\$100 <input type="checkbox"/>	R807	\$100 <input type="checkbox"/>
Butterfly harness	CHOOSE	Regular	R889	\$155 <input type="checkbox"/>	R889	\$155 <input type="checkbox"/>
		Slim cut	R826	\$155 <input type="checkbox"/>	R826	\$155 <input type="checkbox"/>
Thigh belt			R813	\$105 <input type="checkbox"/>	R813	\$105 <input type="checkbox"/>
Pelvic harness			R846	\$200 <input type="checkbox"/>	R846	\$200 <input type="checkbox"/>
Hip guides (pair)			R867	\$221 <input type="checkbox"/>	R867	\$221 <input type="checkbox"/>
Abductor	CHOOSE	Small	R884	\$135 <input type="checkbox"/>	R884	\$135 <input type="checkbox"/>
		Large	R891	\$155 <input type="checkbox"/>	R891	\$155 <input type="checkbox"/>
Adductors (pair)			R892	\$210 <input type="checkbox"/>	R892	\$210 <input type="checkbox"/>
Leg prompt			R859	\$250 <input type="checkbox"/>	R859	\$250 <input type="checkbox"/>
Footboard (recommended for chairs with long legs)			R857	\$520 <input type="checkbox"/>	Included	
Ankle straps (pair)			R855	\$100 <input type="checkbox"/>	R855	\$100 <input type="checkbox"/>
Sandals (pair)			R964	\$115 <input type="checkbox"/>	R964	\$115 <input type="checkbox"/>
Wedges (pair, require sandals)			R699	\$80 <input type="checkbox"/>	R699	\$80 <input type="checkbox"/>
Lumbar and seat support kit			R888	\$95 <input type="checkbox"/>	R888	\$95 <input type="checkbox"/>
Backrest filler pad			R856	\$130 <input type="checkbox"/>	R856	\$130 <input type="checkbox"/>
Whitmyer headrest adapter			R804	\$60 <input type="checkbox"/>	R804	\$60 <input type="checkbox"/>
Total			R860	\$ _____	R870	\$ _____

* Required

Room to grow

Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"



Large chair at its largest dimensions


- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"




Large standard base
R860 Rifton Activity Chair



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R860 large Standard base
Frame width	short legs: 26¾ long legs: 26¾
	short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24
	short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250