

# 2023 order form for medium

Use dimension chart to select appropriate size.



		Standard R840 medium		Hi/lo R850 medium		
* <b>Seat and back</b> includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$640 <input type="checkbox"/>	R851	\$1,335 <input type="checkbox"/>
		Adjustable w/ back spring	R842	\$730 <input type="checkbox"/>	R852	\$1,450 <input type="checkbox"/>
* <b>Base</b>	CHOOSE	Without spring	R843	\$370 <input type="checkbox"/>	R853	\$2,680 <input type="checkbox"/>
		With spring	R844	\$445 <input type="checkbox"/>	N/A	
* <b>Arm supports</b> (pair)	CHOOSE	Armrests	R868	\$260 <input type="checkbox"/>	R868	\$260 <input type="checkbox"/>
		Forearm prompts	R869	\$650 <input type="checkbox"/>	R869	\$650 <input type="checkbox"/>
* <b>Adjustable legs</b>	CHOOSE	Short	R893	\$60 <input type="checkbox"/>	N/A	
		Long	R894	\$85 <input type="checkbox"/>	N/A	
		Short with casters	R895	\$315 <input type="checkbox"/>	N/A	
		Long with casters	R896	\$340 <input type="checkbox"/>	N/A	
* <b>Pads</b>	CHOOSE	Tan	R835	\$260 <input type="checkbox"/>	R835	\$260 <input type="checkbox"/>
		Blue	R836	\$260 <input type="checkbox"/>	R836	\$260 <input type="checkbox"/>
		Red	R837	\$260 <input type="checkbox"/>	R837	\$260 <input type="checkbox"/>
		Pink	R974	\$260 <input type="checkbox"/>	R974	\$260 <input type="checkbox"/>
		Purple	R975	\$260 <input type="checkbox"/>	R975	\$260 <input type="checkbox"/>
		Green	R976	\$260 <input type="checkbox"/>	R976	\$260 <input type="checkbox"/>
<b>Headrest</b>	CHOOSE	Contoured	R897	\$250 <input type="checkbox"/>	R897	\$250 <input type="checkbox"/>
		Flat	R898	\$250 <input type="checkbox"/>	R898	\$250 <input type="checkbox"/>
		Adjustable winged	R881	\$275 <input type="checkbox"/>	R881	\$275 <input type="checkbox"/>
<b>Laterals</b> (pair)	CHOOSE	Small	R828	\$210 <input type="checkbox"/>	R828	\$210 <input type="checkbox"/>
		Small with chest strap	R829	\$260 <input type="checkbox"/>	R829	\$260 <input type="checkbox"/>
		Large	R838	\$225 <input type="checkbox"/>	R838	\$225 <input type="checkbox"/>
<b>Push handles</b> (recommended for all chairs with wheels)			R878	\$160 <input type="checkbox"/>	Included	
<b>Tray</b> (requires armrests)			R668	\$290 <input type="checkbox"/>	R668	\$290 <input type="checkbox"/>
<b>Handhold</b>			R649	\$30 <input type="checkbox"/>	R649	\$30 <input type="checkbox"/>
<b>Chest strap</b> (wide)			R863	\$95 <input type="checkbox"/>	R863	\$95 <input type="checkbox"/>
<b>Butterfly harness</b>	CHOOSE	Regular	R864	\$130 <input type="checkbox"/>	R864	\$130 <input type="checkbox"/>
		Slim cut	R882	\$130 <input type="checkbox"/>	R882	\$130 <input type="checkbox"/>
<b>Thigh belt</b>			R865	\$100 <input type="checkbox"/>	R865	\$100 <input type="checkbox"/>
<b>Pelvic harness</b>			R866	\$180 <input type="checkbox"/>	R866	\$180 <input type="checkbox"/>
<b>Hip guides</b> (pair)			R867	\$221 <input type="checkbox"/>	R867	\$221 <input type="checkbox"/>
<b>Abductor</b>	CHOOSE	Small	R884	\$135 <input type="checkbox"/>	R884	\$135 <input type="checkbox"/>
		Large	R891	\$155 <input type="checkbox"/>	R891	\$155 <input type="checkbox"/>
<b>Adductors</b> (pair)			R892	\$210 <input type="checkbox"/>	R892	\$210 <input type="checkbox"/>
<b>Leg prompt</b>			R858	\$235 <input type="checkbox"/>	R858	\$235 <input type="checkbox"/>
<b>Footboard</b> (recommended for chairs with long legs)			R848	\$505 <input type="checkbox"/>	Included	
<b>Ankle straps</b> (pair)			R899	\$95 <input type="checkbox"/>	R899	\$95 <input type="checkbox"/>
<b>Sandals</b> (pair)			R954	\$100 <input type="checkbox"/>	R954	\$100 <input type="checkbox"/>
<b>Wedges</b> (pair, require sandals)			R698	\$75 <input type="checkbox"/>	R698	\$75 <input type="checkbox"/>
<b>Lumbar and seat support kit</b>			R887	\$85 <input type="checkbox"/>	R887	\$85 <input type="checkbox"/>
<b>Backrest filler pad</b>			R849	\$130 <input type="checkbox"/>	R849	\$130 <input type="checkbox"/>
<b>Whitmyer headrest adapter</b>			R804	\$60 <input type="checkbox"/>	R804	\$60 <input type="checkbox"/>
<b>Total</b>			<b>R840</b>	<b>\$ _____</b>	<b>R850</b>	<b>\$ _____</b>

\* Required

# Room to grow

## Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"



## Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"



**Medium standard base**  
R840 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
<p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	R840 medium Standard base
Frame width	short legs: 23¾
	long legs: 24½
Seat height above floor	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above footboard	short legs: 12½–16½
	long legs: 19–23
Seat angle (tilt-in-space)–degrees	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Backrest angle–degrees	15° forward, 15° back
Footboard knee angle–degrees	5° forward, 20° back
Seat height above footboard	45°–110°
Seat width with hip guides	7¾–15¼
without hip guides	8½–11½
Seat depth	14
Armrest height above seat	11–16
Distance between laterals	7–10½
Backrest height	6½–12
Headrest height above seat	15½–19
Max. working load (lb)	17½–24½
	150



**Medium hi/lo base**  
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
<p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	R850 medium Hi/lo base
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides	8½–11½
without hip guides	14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lb)	150