


2023 order form

Use dimension chart to select appropriate size.



| User dimensions (inches) | E930 small | E940 med. | E950 large |
|---|----------------------------|------------|------------|
| Height | 25 – 48 | 44 – 59 | 57 – 72 |
|  Key user dimension: height Select the appropriate Prone Stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load. | | | |
| | Item dimensions (inches) | E930 small | E940 med. |
| Board length | 23½ (16 without kneeboard) | 30 | 39 |
| Board width (top / bottom) | 12½ / 12½ | 11 / 16 | 14 / 19 |
| Top of board to footboard | 20 – 35 | 33 – 43½ | 44 – 57 |
| Height when horizontal | 29½ | 29½ | 31½ |
| Height when vertical | 38½ | 45½ | 61 |
| Distance between laterals | 8 – 16 | 8 – 16 | 12 – 22 |
| Base width | 23 | 26 | 30 |
| Base length | 35 | 39½ | 44 |
| Tray inside width (widest point) | 22½ | 22½ | 29½ |
| Tray inside length | 17½ | 17½ | 17½ |
| Max. working load (lb) | 100 | 150 | 200 |

* **Prone Stander** includes a pair of trunk laterals, hip strap, padding and footboard.

Specify color: red or blue

| | | E930 small | E940 medium | E950 large |
|---|-------------------------------------|-------------------------------------|---------------------------------------|---------------------------------------|
| CHOOSE | E931 | \$2,300 <input type="checkbox"/> | E941 \$2,670 <input type="checkbox"/> | E951 \$3,270 <input type="checkbox"/> |
| | Red | <input type="checkbox"/> | Red <input type="checkbox"/> | Red <input type="checkbox"/> |
| | Blue | <input type="checkbox"/> | Blue <input type="checkbox"/> | Blue <input type="checkbox"/> |
| Additional trunk laterals (pair) | N/A | E614 \$230 <input type="checkbox"/> | E614 \$230 <input type="checkbox"/> | |
| Knee laterals (pair) | E942 \$215 <input type="checkbox"/> | E942 \$215 <input type="checkbox"/> | E942 \$215 <input type="checkbox"/> | |
| Additional hip strap | N/A | E946 \$85 <input type="checkbox"/> | E956 \$95 <input type="checkbox"/> | |
| Tray | E947 \$645 <input type="checkbox"/> | E947 \$645 <input type="checkbox"/> | E957 \$735 <input type="checkbox"/> | |
| Round abductor | E933 \$110 <input type="checkbox"/> | E933 \$110 <input type="checkbox"/> | E933 \$110 <input type="checkbox"/> | |
| Collar (requires round abductor) | E934 \$90 <input type="checkbox"/> | E934 \$90 <input type="checkbox"/> | E934 \$90 <input type="checkbox"/> | |
| Round abductor with hip stabilizer | E655 \$315 <input type="checkbox"/> | E655 \$315 <input type="checkbox"/> | E655 \$315 <input type="checkbox"/> | |
| Adjustable abduction wedge | N/A | E945 \$260 <input type="checkbox"/> | E945 \$260 <input type="checkbox"/> | |
| Adjustable abduction wedge with hip stabilizer | N/A | E605 \$330 <input type="checkbox"/> | E605 \$330 <input type="checkbox"/> | |
| Sandals (pair) | Small | E422 \$100 <input type="checkbox"/> | E422 \$100 <input type="checkbox"/> | N/A |
| | Medium | E433 \$110 <input type="checkbox"/> | E433 \$110 <input type="checkbox"/> | E433 \$110 <input type="checkbox"/> |
| | Large | N/A | N/A | E432 \$110 <input type="checkbox"/> |
| Wedges (pair, require sandals) | Small | E424 \$70 <input type="checkbox"/> | E424 \$70 <input type="checkbox"/> | N/A |
| | Medium | E435 \$75 <input type="checkbox"/> | E435 \$75 <input type="checkbox"/> | E435 \$75 <input type="checkbox"/> |
| | Large | N/A | N/A | E434 \$80 <input type="checkbox"/> |
| Hand Anchor | K821 \$165 <input type="checkbox"/> | K821 \$165 <input type="checkbox"/> | K821 \$165 <input type="checkbox"/> | |
| Total | E930 \$ _____ | E940 \$ _____ | E950 \$ _____ | |

* Required