

2023 order form for small

Use dimension chart to select appropriate size.



| | | | Standard R820 small | | Hi/lo R830 small | |
|---|--------|---------------------------|------------------------|--------------------------------|---------------------|----------------------------------|
| * Seat and back includes seatbelt (Hi/lo includes footboard and push handles) | CHOOSE | Adjustable | R821 | \$625 <input type="checkbox"/> | R831 | \$1,275 <input type="checkbox"/> |
| | | Adjustable w/ back spring | R822 | \$715 <input type="checkbox"/> | R832 | \$1,390 <input type="checkbox"/> |
| * Base | CHOOSE | Without spring | R823 | \$365 <input type="checkbox"/> | R833 | \$2,180 <input type="checkbox"/> |
| | | With spring | R824 | \$430 <input type="checkbox"/> | N/A | |
| * Arm supports (pair) | CHOOSE | Armrests | R827 | \$255 <input type="checkbox"/> | R827 | \$255 <input type="checkbox"/> |
| | | Forearm prompts | R809 | \$625 <input type="checkbox"/> | R809 | \$625 <input type="checkbox"/> |
| * Adjustable legs | CHOOSE | Short | R803 | \$50 <input type="checkbox"/> | N/A | |
| | | Long | R894 | \$85 <input type="checkbox"/> | N/A | |
| | | Short with casters | R805 | \$310 <input type="checkbox"/> | N/A | |
| | | Long with casters | R896 | \$340 <input type="checkbox"/> | N/A | |
| * Pads | CHOOSE | Tan | R815 | \$255 <input type="checkbox"/> | R815 | \$255 <input type="checkbox"/> |
| | | Blue | R816 | \$255 <input type="checkbox"/> | R816 | \$255 <input type="checkbox"/> |
| | | Red | R817 | \$255 <input type="checkbox"/> | R817 | \$255 <input type="checkbox"/> |
| | | Pink | R971 | \$255 <input type="checkbox"/> | R971 | \$255 <input type="checkbox"/> |
| | | Purple | R972 | \$255 <input type="checkbox"/> | R972 | \$255 <input type="checkbox"/> |
| | | Green | R973 | \$255 <input type="checkbox"/> | R973 | \$255 <input type="checkbox"/> |
| Headrest | CHOOSE | Contoured | R897 | \$250 <input type="checkbox"/> | R897 | \$250 <input type="checkbox"/> |
| | | Flat | R898 | \$250 <input type="checkbox"/> | R898 | \$250 <input type="checkbox"/> |
| | | Adjustable winged | R881 | \$275 <input type="checkbox"/> | R881 | \$275 <input type="checkbox"/> |
| Laterals (pair) | CHOOSE | Small | R828 | \$210 <input type="checkbox"/> | R828 | \$210 <input type="checkbox"/> |
| | | Small with chest strap | R819 | \$255 <input type="checkbox"/> | R819 | \$255 <input type="checkbox"/> |
| Push handles (recommended for all chairs with wheels) | | | R878 | \$160 <input type="checkbox"/> | Included | |
| Tray (requires armrests) | | | R667 | \$255 <input type="checkbox"/> | R667 | \$255 <input type="checkbox"/> |
| Handhold | | | R649 | \$30 <input type="checkbox"/> | R649 | \$30 <input type="checkbox"/> |
| Chest strap (wide) | | | R802 | \$90 <input type="checkbox"/> | R802 | \$90 <input type="checkbox"/> |
| Butterfly harness | | | R834 | \$125 <input type="checkbox"/> | R834 | \$125 <input type="checkbox"/> |
| Thigh belt | | | R811 | \$95 <input type="checkbox"/> | R811 | \$95 <input type="checkbox"/> |
| Pelvic harness | | | R845 | \$160 <input type="checkbox"/> | R845 | \$160 <input type="checkbox"/> |
| Hip guides (pair) | | | R814 | \$215 <input type="checkbox"/> | R814 | \$215 <input type="checkbox"/> |
| Abductor | | | R884 | \$135 <input type="checkbox"/> | R884 | \$135 <input type="checkbox"/> |
| Adductors (pair) | | | R892 | \$210 <input type="checkbox"/> | R892 | \$210 <input type="checkbox"/> |
| Leg prompt | | | R854 | \$225 <input type="checkbox"/> | R854 | \$225 <input type="checkbox"/> |
| Footboard (recommended for chairs with long legs) | | | R847 | \$495 <input type="checkbox"/> | Included | |
| Mini kit (backrest insert and trunk support) | | | R812 | \$135 <input type="checkbox"/> | R812 | \$135 <input type="checkbox"/> |
| Ankle straps (pair) | | | R899 | \$95 <input type="checkbox"/> | R899 | \$95 <input type="checkbox"/> |
| Sandals (pair) | | | R944 | \$95 <input type="checkbox"/> | R944 | \$95 <input type="checkbox"/> |
| Wedges (pair, require sandals) | | | R697 | \$70 <input type="checkbox"/> | R697 | \$70 <input type="checkbox"/> |
| Lumbar and seat support kit | | | R886 | \$70 <input type="checkbox"/> | R886 | \$70 <input type="checkbox"/> |
| Backrest filler pad | | | R849 | \$130 <input type="checkbox"/> | R849 | \$130 <input type="checkbox"/> |
| Whitmyer headrest adapter | | | R804 | \$60 <input type="checkbox"/> | R804 | \$60 <input type="checkbox"/> |
| Total | | | R820 | \$ _____ | R830 | \$ _____ |

* Required

Room to grow

Small chair with mini kit at its smallest dimensions

- 1 Backrest height 12½"
- 2 Trunk support width 5½"
- 3 Seat depth 7"
- 4 Seat width 7"
- 5 Armrest height above seat 5"
- 6 Seat height above footboard 7"




Small chair at its largest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 11½"
- 3 Seat depth 12"
- 4 Seat width 12"
- 5 Armrest height above seat 7½"
- 6 Seat height above footboard 12"




Small standard base
R820 Rifton Activity Chair

| User dimensions (inches) | |
|---|---|
| Height with mini kit | 32–48 28–38 |
|  Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load. | |
| | |
| Item dimensions (inches) | |
| Frame width | R820 small Standard base |
| | short legs: 21 long legs: 23 |
| Seat height above floor | short legs w/ casters: 22½ long legs w/ casters: 23½ |
| | short legs: 9½–12½ long legs: 18½–21½ |
| Seat height above footboard | short legs w/ casters: 13½–16½ long legs w/ casters: 18½–21½ |
| | 15° forward, 15° back |
| Backrest angle—degrees | 5° forward, 20° back |
| Footboard knee angle—degrees | 45°–110° |
| Seat height above footboard | 7–12 |
| Seat width with hip guides | 7–9 |
| without hip guides | 12 |
| Seat depth | 8–12 |
| with mini kit | 7–11 |
| Armrest height above seat | 5–7½ |
| Trunk support width | 5½–11½ |
| Backrest height | 12½–15½ |
| Headrest height above seat | 14½–21 |
| Max. working load (lb) | 75 |



Small hi/lo base
R830 Rifton Activity Chair

| User dimensions (inches) | |
|---|-----------------------|
| Height with mini kit | 32–48 28–38 |
|  Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load. | |
| | |
| Item dimensions (inches) | |
| Frame width | R830 small Hi/lo base |
| | 26 |
| Seat height above floor | 10–25 |
| Seat angle (tilt-in-space)—degrees | 15° forward, 25° back |
| Backrest angle—degrees | 5° forward, 20° back |
| Footboard knee angle—degrees | 45°–110° |
| Seat height above footboard | 7–12 |
| Seat width with hip guides | 7–9 |
| without hip guides | 12 |
| Seat depth | 8–12 |
| with mini kit | 7–11 |
| Armrest height above seat | 5–7½ |
| Trunk support width | 5½–11½ |
| Backrest height | 12½–15½ |
| Headrest height above seat | 14½–21 |
| Max. working load (lb) | 75 |