


How to buy

- 1. Measure** your client to determine the size you need.
(See dimension chart for help.)
- 2. Select** frame and color.
- 3. Choose** the accessories you need.
- 4. Add** up the total cost.
- 5. Order** online in minutes at **Rifton.com**.
- 6. Call** us any time at **800.571.8198** if you need assistance.
We're always here to help.

User dimensions (inches)	K509 mini	K501 small	
Elbow height	15½ – 20½	18½ – 27½	
 <p>Key user dimension: elbow height Measure the vertical distance from bent elbow to floor while the user is standing upright. Choose the gait trainer that allows for growth.</p> <p>Important: Make sure that chest prompt width allows for growth. User's weight must not exceed the maximum working load.</p> <p>Small or medium prompts are recommended for most users of the medium Pacer.</p>			
	Item dimensions (inches)	K509 mini	K501 small
	Arm prompt height	15½ – 20½	18½ – 27½
Overall width	20½	22½	
Overall length	22½	27	
Chest prompt height (top edge)	17½ – 22½	22 – 30½	
Frame height	12½	15½ – 19½	
Frame weight (lbs)	7½	11½	
Max. working load (lbs)	50	75	



Pacer colors



Red



Blue



Lime



Raspberry

	K509 mini		K501 small	
Frame	K510	\$655 <input type="checkbox"/>	K511	\$885 <input type="checkbox"/>
Specify color: Red, blue, lime or raspberry	CHOOSE	RED <input type="checkbox"/>	RAS <input type="checkbox"/>	RED <input type="checkbox"/>
		BLU <input type="checkbox"/>	LIM <input type="checkbox"/>	BLU <input type="checkbox"/>
				LIM <input type="checkbox"/>

		K509 mini		K501 small	
Optional accessories	Arm prompts (pair)	K512	\$440 <input type="checkbox"/>	K512	\$440 <input type="checkbox"/>
	Hand loops	CHOOSE	Single	K578 \$88 <input type="checkbox"/>	K578 \$88 <input type="checkbox"/>
			Pair	K579 \$170 <input type="checkbox"/>	K579 \$170 <input type="checkbox"/>
	Chest prompt	CHOOSE	Small	K514 \$490 <input type="checkbox"/>	K514 \$490 <input type="checkbox"/>
			Medium	N/A	K524 \$530 <input type="checkbox"/>
	Pelvic support (includes handholds)	CHOOSE	Small	K556 \$360 <input type="checkbox"/>	K556 \$360 <input type="checkbox"/>
			Medium	N/A	K557 \$380 <input type="checkbox"/>
	Hip positioner (includes handholds)	K513	\$360 <input type="checkbox"/>	K513	\$360 <input type="checkbox"/>
	Hip positioner pad	K541	\$95 <input type="checkbox"/>	K541	\$95 <input type="checkbox"/>
	Thigh prompts (pair)	N/A		K515	\$365 <input type="checkbox"/>
	Ankle prompts (pair)	K529	\$145 <input type="checkbox"/>	K516	\$155 <input type="checkbox"/>
	Additional handholds (pair)	K519	\$98 <input type="checkbox"/>	K519	\$98 <input type="checkbox"/>
	Communication tray	K518	\$335 <input type="checkbox"/>	K518	\$335 <input type="checkbox"/>
Attendant guide bar	K589	\$125 <input type="checkbox"/>	K589	\$125 <input type="checkbox"/>	
Accessories tote	K522	\$85 <input type="checkbox"/>	K522	\$85 <input type="checkbox"/>	
Total	K509	\$ _____	K501	\$ _____	