

Pacer Positioning Checklist

Client Name: _____

Therapist Name: _____

Date: _____

Use this Positioning Checklist as a convenient way to assure optimal use of the Pacer. Jot down notes to customize your instructions for each individual.

Frame Height



Top bar of frame at level of:

- Hip trochanters
- Waist level
- Above waist


Notes: _____

Chest Prompt



- Clamp Position: Inside frame Outside frame
- Width Adjustment: Snug to Trunk OR Width measurement: _____
- Height Adjustment: Specify # of notches visible above clamp: _____

Left Arm Prompt




- Clamp Placement: Front Bar Side Bar
- Clamp Position: Inside frame Outside frame
- L-bar Support: Points Forward Points Backward
 Points toward R Points toward L
- L-bar Height: Specify # of notches visible above clamp: _____
- Arm Prompt: Angle Up _____ Angle Down _____
Rotation Position (specify as o'clock when standing in Pacer) _____

Right Arm Prompt



- Clamp Placement: Front Bar Side Bar
- Clamp Position: Inside frame Outside frame
- L-bar Support: Points Forward Points Backward
 Points toward R Points toward L
- L-bar Height: Specify # of notches visible above clamp: _____
- Arm Prompt: Angle Up _____ Angle Down _____
Rotation Position (specify as o'clock when standing in Pacer) _____

Hip Positioner



2 Front Straps: Explain location where attached: _____

Handhold Height (for 2 Back Straps): Specify # of notches visible above clamp: _____

Handhold Clamps: Specify distance from end of top bar: _____

Notes: _____

Thigh Prompts



Right: Clamp Inside Clamp Outside

Left: Clamp Inside Clamp Outside

Explain Position: Right _____

Left _____

Ankle Prompts



Describe location of Spring Adjusters: Right _____

Left _____

Casters



Front Casters: Set Parallel Set Ratchet Set Drag

Rear Casters: Set Parallel Set Ratchet Set Drag