

# More comfort with every step.

**The Pacer's standard, rigid hip positioner isn't for everyone.** That's why we developed a narrower, less restrictive alternative – the **pelvic support**.

**Like the hip positioner**, it offers temporary weight-bearing assistance as needed, while the user builds strength and stamina. But the **pelvic support has a softer design** that provides a **more comfortable fit** for larger users, those whose legs have a tendency to scissor, or anyone for whom the hip positioner is simply too bulky.

- **Less chafing** for larger users and those whose gait makes the hip positioner uncomfortable
- **Narrower design**, ideal for the youngest walkers and others who need minimal abduction
- **Soft and strong** foam and web-reinforced construction is comfortable, yet durable



“ The new pelvic support for the Pacer has worked amazingly well for a student who had high post-operative discomfort from hip de-rotation and foot osteotomies. He had been unable to do much weight-bearing, but with the new pelvic support, he doubled his walking distance on the first day! Thank you! ”

Jennifer Lynch DPT  
Monroe-Woodbury  
Central School District, NY

New pelvic supports  
for Rifton Pacers



Products for life from people who care.™

# Simply choose the size pelvic support to fit your user and Pacer.

## Three sizes of Pelvic Supports

Pacer Sizes >>

Mini

Small

Medium

Large

XL



### Small

(Max. working load 75 lbs)  
without handholds K553 **\$245**  
(see step 2 below)  
with handholds K556 **\$335**



### Medium

(Max. working load 275 lbs)  
without handholds K554 **\$265**  
(see step 2 below)  
with handholds K557 **\$355**



### Large

(Max. working load 275 lbs)  
without handholds K555 **\$280**  
(see step 2 below)  
with handholds K558 **\$370**



### Key user dimension:

**One-half the thigh circumference**

Small pelvic support: 4" – 8"

Medium pelvic support: 6" – 10"

Large pelvic support: 8" – 14"



The pelvic support is made of microbe-resistant foam padding sandwiched between two layers of durable polyurethane. Routine cleaning is easy with a 10% bleach solution.



### Choose the **hip positioner** for:

- More abduction
- More rigid positioning

### Choose the **pelvic support** for:

- Minimal abduction
- More freedom of movement

## How to buy

**1. Measure** your client's thigh circumference and divide in half to determine the size of pelvic support you need. (Refer to the charts above.)

**2. Determine** if you already have two handholds (with clamps) to dedicate to the pelvic support – if not, be sure to include them with the order.

**3. Add** up the total cost.

**4. Order** online in minutes at **Rifton.com**. Or call us at **800.571.8198** to order or for assistance.



For some people, our standard hip positioner is too restrictive or provides too much abduction. For these clients, the pelvic support will provide all the positioning and weight bearing needed, with none of the bulk and rigidity of our standard hip positioner.

Call 800.571.8198 or go to **Rifton.com**



Products for life from people who care.™