


Compass Chair dimensions

User dimensions (inches)	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
Height	36 – 45	43 – 50	48 – 57	55 – 65	65+
Knee to heel	With short legs: 9 – 11 With standard legs: 10 – 13	10½ – 14½	12½ – 16½	14½ – 18½	16½ – 20½
Min – max thigh circumference for	10½ – 23½	12 – 28½	14½ – 30½	15 – 34½	17½ – 38½
 <p>Key user dimension: height / knee to heel The user's overall height and measurement from knee to heel is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each user, and that the user's weight does not exceed the maximum working load.</p>					
Item dimensions (inches)	R310	R320	R330	R340	R350
Seat height – in 1" increments	9 – 13	10½ – 14½	12½ – 16½	14½ – 18½	16½ – 20½
Seat depth	9½	11	12¾	14¾	16½
Seat width	10½	12¾	14	16¾	18¾
Backrest height	9½	11	12½	14¾	16
Armrest height above the seat	5½	6	7	8¾	8¾
Armrest length	6½	7½	8½	9¾	11
Max. working load (lb)	100	120	175	220	250

2025 order form

Scan for sizing and pricing:
<https://www.rifton.com/compasschair>



Use dimension chart to select appropriate size.

		R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
* Compass chair		R311 \$200 <input type="checkbox"/>	R321 \$230 <input type="checkbox"/>	R331 \$260 <input type="checkbox"/>	R341 \$290 <input type="checkbox"/>	R351 \$320 <input type="checkbox"/>
* Leg type	Short	R314 \$155 <input type="checkbox"/>	N/A	N/A	N/A	N/A
	Short (low friction)	R334 \$155 <input type="checkbox"/>	N/A	N/A	N/A	N/A
	Standard	R315* \$155 <input type="checkbox"/>	R315 \$155 <input type="checkbox"/>	R315 \$155 <input type="checkbox"/>	R315 \$155 <input type="checkbox"/>	R315 \$155 <input type="checkbox"/>
	Standard (low friction)	R335* \$155 <input type="checkbox"/>	R335 \$155 <input type="checkbox"/>	R335 \$155 <input type="checkbox"/>	R335 \$155 <input type="checkbox"/>	R335 \$155 <input type="checkbox"/>
New! Front stability feet (pair)		R326 \$40 <input type="checkbox"/>	R326 \$40 <input type="checkbox"/>	R326 \$40 <input type="checkbox"/>	R326 \$40 <input type="checkbox"/>	R326 \$40 <input type="checkbox"/>
New! Rear stability feet (pair)		R327 \$40 <input type="checkbox"/>	R327 \$40 <input type="checkbox"/>	R327 \$40 <input type="checkbox"/>	R327 \$40 <input type="checkbox"/>	R327 \$40 <input type="checkbox"/>
New! Dynamic feet (set of four)		R328 \$100 <input type="checkbox"/>	R328 \$100 <input type="checkbox"/>	R328 \$100 <input type="checkbox"/>	R328 \$100 <input type="checkbox"/>	R328 \$100 <input type="checkbox"/>
Seat belt		R312 \$48 <input type="checkbox"/>	R312 \$48 <input type="checkbox"/>	R332 \$48 <input type="checkbox"/>	R342 \$48 <input type="checkbox"/>	R342 \$48 <input type="checkbox"/>
Pelvic harness		R313 \$180 <input type="checkbox"/>	R323 \$190 <input type="checkbox"/>	R333 \$200 <input type="checkbox"/>	R343 \$210 <input type="checkbox"/>	R353 \$220 <input type="checkbox"/>

*** Required**

* The standard legs can be used on size 1, but the lower two positions will not be accessible.