

# 2026 order form for large

Use dimension chart to select appropriate size.

Scan for sizing and pricing:  
<https://www.rifton.com/activitychair>



		Standard R860 large			Hi/lo R870 large			
* <b>Seat and back</b> includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861	\$820	<input type="checkbox"/>	R871	\$1,550	<input type="checkbox"/>
		Adjustable w/ back spring	R862	\$940	<input type="checkbox"/>	R872	\$1,680	<input type="checkbox"/>
* <b>Base</b>	CHOOSE	Without spring	R873	\$495	<input type="checkbox"/>	R883	\$3,470	<input type="checkbox"/>
		With spring	R874	\$600	<input type="checkbox"/>	N/A		
* <b>Arm supports</b> (pair)	CHOOSE	Armrests	R868	\$290	<input type="checkbox"/>	R868	\$290	<input type="checkbox"/>
		Forearm prompts	R879	\$750	<input type="checkbox"/>	R879	\$750	<input type="checkbox"/>
* <b>Adjustable legs</b>	CHOOSE	Short	R808	\$83	<input type="checkbox"/>	N/A		
		Long	R894	\$93	<input type="checkbox"/>	N/A		
		Short with casters	R895	\$355	<input type="checkbox"/>	N/A		
		Long with casters	R896	\$390	<input type="checkbox"/>	N/A		
* <b>Pads</b>	CHOOSE	Tan	R875	\$310	<input type="checkbox"/>	R875	\$310	<input type="checkbox"/>
		Blue	R876	\$310	<input type="checkbox"/>	R876	\$310	<input type="checkbox"/>
		Red	R877	\$310	<input type="checkbox"/>	R877	\$310	<input type="checkbox"/>
		Pink	R977	\$310	<input type="checkbox"/>	R977	\$310	<input type="checkbox"/>
		Purple	R978	\$310	<input type="checkbox"/>	R978	\$310	<input type="checkbox"/>
		Green	R979	\$310	<input type="checkbox"/>	R979	\$310	<input type="checkbox"/>
<b>Headrest</b>	CHOOSE	Contoured	R897	\$280	<input type="checkbox"/>	R897	\$280	<input type="checkbox"/>
		Flat	R898	\$280	<input type="checkbox"/>	R898	\$280	<input type="checkbox"/>
		Adjustable winged	R881	\$310	<input type="checkbox"/>	R881	\$310	<input type="checkbox"/>
<b>Laterals</b> (pair)	CHOOSE	Small	R828	\$240	<input type="checkbox"/>	R828	\$240	<input type="checkbox"/>
		Large	R838	\$255	<input type="checkbox"/>	R838	\$255	<input type="checkbox"/>
		Large with chest strap	R839	\$305	<input type="checkbox"/>	R839	\$305	<input type="checkbox"/>
<b>Push handles</b> (recommended for all chairs with wheels)		R878	\$175	<input type="checkbox"/>	Included			
<b>Tray</b> (requires armrests)		R669	\$345	<input type="checkbox"/>	R669	\$345	<input type="checkbox"/>	
<b>Handhold</b>		R649	\$33	<input type="checkbox"/>	R649	\$33	<input type="checkbox"/>	
<b>Chest strap</b> (wide)		R807	\$115	<input type="checkbox"/>	R807	\$115	<input type="checkbox"/>	
<b>Butterfly harness</b>	CHOOSE	Regular	R889	\$175	<input type="checkbox"/>	R889	\$175	<input type="checkbox"/>
		Slim cut	R826	\$175	<input type="checkbox"/>	R826	\$175	<input type="checkbox"/>
<b>Thigh belt</b>		R813	\$120	<input type="checkbox"/>	R813	\$120	<input type="checkbox"/>	
<b>Pelvic harness</b>		R846	\$230	<input type="checkbox"/>	R846	\$230	<input type="checkbox"/>	
<b>Hip guides</b> (pair)		R867	\$250	<input type="checkbox"/>	R867	\$250	<input type="checkbox"/>	
<b>Abductor</b>	CHOOSE	Small	R884	\$150	<input type="checkbox"/>	R884	\$150	<input type="checkbox"/>
		Large	R891	\$170	<input type="checkbox"/>	R891	\$170	<input type="checkbox"/>
<b>Adductors</b> (pair)		R892	\$240	<input type="checkbox"/>	R892	\$240	<input type="checkbox"/>	
<b>Leg prompt</b>		R859	\$280	<input type="checkbox"/>	R859	\$280	<input type="checkbox"/>	
<b>Footboard</b> (recommended for chairs with long legs)		R857	\$580	<input type="checkbox"/>	Included			
<b>Ankle straps</b> (pair)		R855	\$115	<input type="checkbox"/>	R855	\$115	<input type="checkbox"/>	
<b>Sandals</b> (pair)		R964	\$130	<input type="checkbox"/>	R964	\$130	<input type="checkbox"/>	
<b>Wedges</b> (pair, require sandals)		R699	\$95	<input type="checkbox"/>	R699	\$95	<input type="checkbox"/>	
<b>Lumbar and seat support kit</b>		R888	\$110	<input type="checkbox"/>	R888	\$110	<input type="checkbox"/>	
<b>Backrest filler pad</b>		R856	\$145	<input type="checkbox"/>	R856	\$145	<input type="checkbox"/>	
<b>Whitmyer headrest adapter</b>		R804	\$67	<input type="checkbox"/>	R804	\$67	<input type="checkbox"/>	
<b>Universal device mount</b>		R967	\$135	<input type="checkbox"/>	R967	\$135	<input type="checkbox"/>	
<b>Directional assist wheel retrofit kit</b>		N/A			R965	\$300	<input type="checkbox"/>	R966 (5 kits) \$1,500 <input type="checkbox"/>

\* Required



# Room to grow

## Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"



## Large chair at its largest dimensions


- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"




**Large standard base**  
R860 Rifton Activity Chair



**Large hi/lo base**  
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R860 large Standard base
Frame width	short legs: 26¾ long legs: 26¾ short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24 short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250

User dimensions (inches)	
Height	50–74
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250