

# 2025 order form for large

Use dimension chart to select appropriate size.

Scan for sizing and pricing:  
<https://www.rifton.com/activitychair>



		Standard R860 large	Hi/lo R870 large
* <b>Seat and back</b> includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861 \$790 <input type="checkbox"/> R871 \$1,500 <input type="checkbox"/>
		Adjustable w/ back spring	R862 \$910 <input type="checkbox"/> R872 \$1,625 <input type="checkbox"/>
* <b>Base</b>	CHOOSE	Without spring	R873 \$480 <input type="checkbox"/> R883 \$3,350 <input type="checkbox"/>
		With spring	R874 \$580 <input type="checkbox"/> N/A
* <b>Arm supports</b> (pair)	CHOOSE	Armrests	R868 \$280 <input type="checkbox"/> R868 \$280 <input type="checkbox"/>
		Forearm prompts	R879 \$725 <input type="checkbox"/> R879 \$725 <input type="checkbox"/>
* <b>Adjustable legs</b>	CHOOSE	Short	R808 \$80 <input type="checkbox"/> N/A
		Long	R894 \$90 <input type="checkbox"/> N/A
		Short with casters	R895 \$340 <input type="checkbox"/> N/A
		Long with casters	R896 \$375 <input type="checkbox"/> N/A
* <b>Pads</b>	CHOOSE	Tan	R875 \$295 <input type="checkbox"/> R875 \$295 <input type="checkbox"/>
		Blue	R876 \$295 <input type="checkbox"/> R876 \$295 <input type="checkbox"/>
		Red	R877 \$295 <input type="checkbox"/> R877 \$295 <input type="checkbox"/>
		Pink	R977 \$295 <input type="checkbox"/> R977 \$295 <input type="checkbox"/>
		Purple	R978 \$295 <input type="checkbox"/> R978 \$295 <input type="checkbox"/>
		Green	R979 \$295 <input type="checkbox"/> R979 \$295 <input type="checkbox"/>
<b>Headrest</b>	CHOOSE	Contoured	R897 \$270 <input type="checkbox"/> R897 \$270 <input type="checkbox"/>
		Flat	R898 \$270 <input type="checkbox"/> R898 \$270 <input type="checkbox"/>
		Adjustable winged	R881 \$300 <input type="checkbox"/> R881 \$300 <input type="checkbox"/>
<b>Laterals</b> (pair)	CHOOSE	Small	R828 \$230 <input type="checkbox"/> R828 \$230 <input type="checkbox"/>
		Large	R838 \$244 <input type="checkbox"/> R838 \$244 <input type="checkbox"/>
		Large with chest strap	R839 \$295 <input type="checkbox"/> R839 \$295 <input type="checkbox"/>
<b>Push handles</b> (recommended for all chairs with wheels)		R878 \$170 <input type="checkbox"/>	Included
<b>Tray</b> (requires armrests)		R669 \$330 <input type="checkbox"/>	R669 \$330 <input type="checkbox"/>
<b>Handhold</b>		R649 \$32 <input type="checkbox"/>	R649 \$32 <input type="checkbox"/>
<b>Chest strap</b> (wide)		R807 \$110 <input type="checkbox"/>	R807 \$110 <input type="checkbox"/>
<b>Butterfly harness</b>	CHOOSE	Regular	R889 \$170 <input type="checkbox"/> R889 \$170 <input type="checkbox"/>
		Slim cut	R826 \$170 <input type="checkbox"/> R826 \$170 <input type="checkbox"/>
<b>Thigh belt</b>		R813 \$115 <input type="checkbox"/>	R813 \$115 <input type="checkbox"/>
<b>Pelvic harness</b>		R846 \$220 <input type="checkbox"/>	R846 \$220 <input type="checkbox"/>
<b>Hip guides</b> (pair)		R867 \$240 <input type="checkbox"/>	R867 \$240 <input type="checkbox"/>
<b>Abductor</b>	CHOOSE	Small	R884 \$145 <input type="checkbox"/> R884 \$145 <input type="checkbox"/>
		Large	R891 \$165 <input type="checkbox"/> R891 \$165 <input type="checkbox"/>
<b>Adductors</b> (pair)		R892 \$230 <input type="checkbox"/>	R892 \$230 <input type="checkbox"/>
<b>Leg prompt</b>		R859 \$270 <input type="checkbox"/>	R859 \$270 <input type="checkbox"/>
<b>Footboard</b> (recommended for chairs with long legs)		R857 \$560 <input type="checkbox"/>	Included
<b>Ankle straps</b> (pair)		R855 \$110 <input type="checkbox"/>	R855 \$110 <input type="checkbox"/>
<b>Sandals</b> (pair)		R964 \$125 <input type="checkbox"/>	R964 \$125 <input type="checkbox"/>
<b>Wedges</b> (pair, require sandals)		R699 \$90 <input type="checkbox"/>	R699 \$90 <input type="checkbox"/>
<b>Lumbar and seat support kit</b>		R888 \$105 <input type="checkbox"/>	R888 \$105 <input type="checkbox"/>
<b>Backrest filler pad</b>		R856 \$140 <input type="checkbox"/>	R856 \$140 <input type="checkbox"/>
<b>Whitmyer headrest adapter</b>		R804 \$65 <input type="checkbox"/>	R804 \$65 <input type="checkbox"/>
<b>Directional Assist Wheel retrofit kit</b>		N/A	R965 \$290 <input type="checkbox"/> R966 (5 kits) \$1450 <input type="checkbox"/>

\* Required



# Room to grow

## Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"



## Large chair at its largest dimensions

- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"



**Large standard base**  
R860 Rifton Activity Chair



**Large hi/lo base**  
R870 Rifton Activity Chair

### User dimensions (inches)

Height 50–74



**Key user dimension: height**

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

**Important:** Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

### Item dimensions (inches)

#### R860 large Standard base

Frame width	short legs: 26¾ long legs: 26¾ short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24 short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250

### User dimensions (inches)

Height 50–74



**Key user dimension: height**

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

**Important:** Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

### Item dimensions (inches)

#### R870 large Hi/lo base

Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250