

# 2026 order form for medium

Use dimension chart to select appropriate size.

Scan for sizing and pricing:  
<https://www.rifton.com/activitychair>



		Standard R840 medium			Hi/lo R850 medium			
* <b>Seat and back</b> includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$720	<input type="checkbox"/>	R851	\$1,490	<input type="checkbox"/>
		Adjustable w/ back spring	R842	\$820	<input type="checkbox"/>	R852	\$1,620	<input type="checkbox"/>
* <b>Base</b>	CHOOSE	Without spring	R843	\$415	<input type="checkbox"/>	R853	\$3,050	<input type="checkbox"/>
		With spring	R844	\$495	<input type="checkbox"/>	N/A		
* <b>Arm supports</b> (pair)	CHOOSE	Armrests	R868	\$290	<input type="checkbox"/>	R868	\$290	<input type="checkbox"/>
		Forearm prompts	R869	\$725	<input type="checkbox"/>	R869	\$725	<input type="checkbox"/>
* <b>Adjustable legs</b>	CHOOSE	Short	R893	\$72	<input type="checkbox"/>	N/A		
		Long	R894	\$93	<input type="checkbox"/>	N/A		
		Short with casters	R895	\$355	<input type="checkbox"/>	N/A		
		Long with casters	R896	\$390	<input type="checkbox"/>	N/A		
* <b>Pads</b>	CHOOSE	Tan	R835	\$300	<input type="checkbox"/>	R835	\$300	<input type="checkbox"/>
		Blue	R836	\$300	<input type="checkbox"/>	R836	\$300	<input type="checkbox"/>
		Red	R837	\$300	<input type="checkbox"/>	R837	\$300	<input type="checkbox"/>
		Pink	R974	\$300	<input type="checkbox"/>	R974	\$300	<input type="checkbox"/>
		Purple	R975	\$300	<input type="checkbox"/>	R975	\$300	<input type="checkbox"/>
		Green	R976	\$300	<input type="checkbox"/>	R976	\$300	<input type="checkbox"/>
<b>Headrest</b>	CHOOSE	Contoured	R897	\$280	<input type="checkbox"/>	R897	\$280	<input type="checkbox"/>
		Flat	R898	\$280	<input type="checkbox"/>	R898	\$280	<input type="checkbox"/>
		Adjustable winged	R881	\$310	<input type="checkbox"/>	R881	\$310	<input type="checkbox"/>
<b>Laterals</b> (pair)	CHOOSE	Small	R828	\$240	<input type="checkbox"/>	R828	\$240	<input type="checkbox"/>
		Small with chest strap	R829	\$295	<input type="checkbox"/>	R829	\$295	<input type="checkbox"/>
		Large	R838	\$255	<input type="checkbox"/>	R838	\$255	<input type="checkbox"/>
<b>Push handles</b> (recommended for all chairs with wheels)		R878	\$175	<input type="checkbox"/>	Included			
<b>Tray</b> (requires armrests)		R668	\$320	<input type="checkbox"/>	R668	\$320	<input type="checkbox"/>	
<b>Handhold</b>		R649	\$33	<input type="checkbox"/>	R649	\$33	<input type="checkbox"/>	
<b>Chest strap</b> (wide)		R863	\$110	<input type="checkbox"/>	R863	\$110	<input type="checkbox"/>	
<b>Butterfly harness</b>	CHOOSE	Regular	R864	\$145	<input type="checkbox"/>	R864	\$145	<input type="checkbox"/>
		Slim cut	R882	\$145	<input type="checkbox"/>	R882	\$145	<input type="checkbox"/>
<b>Thigh belt</b>		R865	\$115	<input type="checkbox"/>	R865	\$115	<input type="checkbox"/>	
<b>Pelvic harness</b>		R866	\$210	<input type="checkbox"/>	R866	\$210	<input type="checkbox"/>	
<b>Hip guides</b> (pair)		R867	\$250	<input type="checkbox"/>	R867	\$250	<input type="checkbox"/>	
<b>Abductor</b>	CHOOSE	Small	R884	\$150	<input type="checkbox"/>	R884	\$150	<input type="checkbox"/>
		Large	R891	\$170	<input type="checkbox"/>	R891	\$170	<input type="checkbox"/>
<b>Adductors</b> (pair)		R892	\$240	<input type="checkbox"/>	R892	\$240	<input type="checkbox"/>	
<b>Leg prompt</b>		R858	\$265	<input type="checkbox"/>	R858	\$265	<input type="checkbox"/>	
<b>Footboard</b> (recommended for chairs with long legs)		R848	\$565	<input type="checkbox"/>	Included			
<b>Ankle straps</b> (pair)		R899	\$110	<input type="checkbox"/>	R899	\$110	<input type="checkbox"/>	
<b>Sandals</b> (pair)		R954	\$115	<input type="checkbox"/>	R954	\$115	<input type="checkbox"/>	
<b>Wedges</b> (pair, require sandals)		R698	\$90	<input type="checkbox"/>	R698	\$90	<input type="checkbox"/>	
<b>Lumbar and seat support kit</b>		R887	\$93	<input type="checkbox"/>	R887	\$93	<input type="checkbox"/>	
<b>Backrest filler pad</b>		R849	\$145	<input type="checkbox"/>	R849	\$145	<input type="checkbox"/>	
<b>Whitmyer headrest adapter</b>		R804	\$67	<input type="checkbox"/>	R804	\$67	<input type="checkbox"/>	
<b>Universal device mount</b>		R967	\$135	<input type="checkbox"/>	R967	\$135	<input type="checkbox"/>	
<b>Directional assist wheel retrofit kit</b>		N/A			R965	\$300 <input type="checkbox"/>	R966 (5 kits) \$1,500 <input type="checkbox"/>	

\* Required

# Room to grow

## Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"



## Medium chair at its largest dimensions


- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




**Medium standard base**  
R840 Rifton Activity Chair



**Medium hi/lo base**  
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R840 medium Standard base
Frame width	short legs: 23¾
	long legs: 24½
Seat height above floor	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above footboard	short legs: 12½–16½
	long legs: 19–23
Seat angle (tilt-in-space)–degrees	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Backrest angle–degrees	15° forward, 15° back
Footboard knee angle–degrees	5° forward, 20° back
Seat height above footboard	45°–110°
Seat width with hip guides without hip guides	7¾–15¼
Seat depth	8½–11½
Armrest height above seat	14
Distance between laterals	11–16
Backrest height	7–10½
Headrest height above seat	6½–12
Max. working load (lb)	15½–19
	17½–24½
	150

User dimensions (inches)	
Height	40–56
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R850 medium Hi/lo base
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides without hip guides	8½–11½
Seat depth	14
Armrest height above seat	11–16
Distance between laterals	7–10½
Backrest height	6½–12
Headrest height above seat	15½–19
Max. working load (lb)	17½–24½
	150