

# 2025 order form for medium

Use dimension chart to select appropriate size.

Scan for sizing and pricing:  
<https://www.rifton.com/activitychair>



		Standard R840 medium			Hi/lo R850 medium						
<b>* Seat and back</b> includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$695	<input type="checkbox"/>	R851	\$1,440	<input type="checkbox"/>			
		Adjustable w/ back spring	R842	\$790	<input type="checkbox"/>	R852	\$1,565	<input type="checkbox"/>			
<b>* Base</b>	CHOOSE	Without spring	R843	\$400	<input type="checkbox"/>	R853	\$2,950	<input type="checkbox"/>			
		With spring	R844	\$480	<input type="checkbox"/>	N/A					
<b>* Arm supports</b> (pair)	CHOOSE	Armrests	R868	\$280	<input type="checkbox"/>	R868	\$280	<input type="checkbox"/>			
		Forearm prompts	R869	\$700	<input type="checkbox"/>	R869	\$700	<input type="checkbox"/>			
<b>* Adjustable legs</b>	CHOOSE	Short	R893	\$70	<input type="checkbox"/>	N/A					
		Long	R894	\$90	<input type="checkbox"/>	N/A					
		Short with casters	R895	\$340	<input type="checkbox"/>	N/A					
		Long with casters	R896	\$375	<input type="checkbox"/>	N/A					
<b>* Pads</b>	CHOOSE	Tan	R835	\$285	<input type="checkbox"/>	R835	\$285	<input type="checkbox"/>			
		Blue	R836	\$285	<input type="checkbox"/>	R836	\$285	<input type="checkbox"/>			
		Red	R837	\$285	<input type="checkbox"/>	R837	\$285	<input type="checkbox"/>			
		Pink	R974	\$285	<input type="checkbox"/>	R974	\$285	<input type="checkbox"/>			
		Purple	R975	\$285	<input type="checkbox"/>	R975	\$285	<input type="checkbox"/>			
		Green	R976	\$285	<input type="checkbox"/>	R976	\$285	<input type="checkbox"/>			
<b>Headrest</b>	CHOOSE	Contoured	R897	\$270	<input type="checkbox"/>	R897	\$270	<input type="checkbox"/>			
		Flat	R898	\$270	<input type="checkbox"/>	R898	\$270	<input type="checkbox"/>			
		Adjustable winged	R881	\$300	<input type="checkbox"/>	R881	\$300	<input type="checkbox"/>			
<b>Laterals</b> (pair)	CHOOSE	Small	R828	\$230	<input type="checkbox"/>	R828	\$230	<input type="checkbox"/>			
		Small with chest strap	R829	\$285	<input type="checkbox"/>	R829	\$285	<input type="checkbox"/>			
		Large	R838	\$244	<input type="checkbox"/>	R838	\$244	<input type="checkbox"/>			
<b>Push handles</b> (recommended for all chairs with wheels)			R878	\$170	<input type="checkbox"/>	Included					
<b>Tray</b> (requires armrests)			R668	\$310	<input type="checkbox"/>	R668	\$310	<input type="checkbox"/>			
<b>Handhold</b>			R649	\$32	<input type="checkbox"/>	R649	\$32	<input type="checkbox"/>			
<b>Chest strap</b> (wide)			R863	\$105	<input type="checkbox"/>	R863	\$105	<input type="checkbox"/>			
<b>Butterfly harness</b>	CHOOSE	Regular	R864	\$140	<input type="checkbox"/>	R864	\$140	<input type="checkbox"/>			
		Slim cut	R882	\$140	<input type="checkbox"/>	R882	\$140	<input type="checkbox"/>			
<b>Thigh belt</b>			R865	\$110	<input type="checkbox"/>	R865	\$110	<input type="checkbox"/>			
<b>Pelvic harness</b>			R866	\$200	<input type="checkbox"/>	R866	\$200	<input type="checkbox"/>			
<b>Hip guides</b> (pair)			R867	\$240	<input type="checkbox"/>	R867	\$240	<input type="checkbox"/>			
<b>Abductor</b>	CHOOSE	Small	R884	\$145	<input type="checkbox"/>	R884	\$145	<input type="checkbox"/>			
		Large	R891	\$165	<input type="checkbox"/>	R891	\$165	<input type="checkbox"/>			
<b>Adductors</b> (pair)			R892	\$230	<input type="checkbox"/>	R892	\$230	<input type="checkbox"/>			
<b>Leg prompt</b>			R858	\$255	<input type="checkbox"/>	R858	\$255	<input type="checkbox"/>			
<b>Footboard</b> (recommended for chairs with long legs)			R848	\$545	<input type="checkbox"/>	Included					
<b>Ankle straps</b> (pair)			R899	\$105	<input type="checkbox"/>	R899	\$105	<input type="checkbox"/>			
<b>Sandals</b> (pair)			R954	\$110	<input type="checkbox"/>	R954	\$110	<input type="checkbox"/>			
<b>Wedges</b> (pair, require sandals)			R698	\$85	<input type="checkbox"/>	R698	\$85	<input type="checkbox"/>			
<b>Lumbar and seat support kit</b>			R887	\$90	<input type="checkbox"/>	R887	\$90	<input type="checkbox"/>			
<b>Backrest filler pad</b>			R849	\$140	<input type="checkbox"/>	R849	\$140	<input type="checkbox"/>			
<b>Whitmyer headrest adapter</b>			R804	\$65	<input type="checkbox"/>	R804	\$65	<input type="checkbox"/>			
<b>Universal device mount</b>			R967	\$130	<input type="checkbox"/>	R967	\$130	<input type="checkbox"/>			
<b>Directional assist wheel retrofit kit</b>			N/A			R965	\$290	<input type="checkbox"/>	R966 (5 kits)	\$1,450	<input type="checkbox"/>

\* Required

# Room to grow

## Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"



## Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




**Medium standard base**  
R840 Rifton Activity Chair



**Medium hi/lo base**  
R850 Rifton Activity Chair

### User dimensions (inches)

Height 40–56



#### Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

**Important:** Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

### Item dimensions (inches)

#### R840 medium Standard base

Frame width	short legs: 23¼ long legs: 24½ short legs w/ casters: 24½ long legs w/ casters: 25½
Seat height above floor	short legs: 12½–16½ long legs: 19–23 short legs w/ casters: 16–20 long legs w/ casters: 19–23
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides without hip guides	8½–11½ 14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lb)	150

### User dimensions (inches)

Height 40–56



#### Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

**Important:** Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

### Item dimensions (inches)

#### R850 medium Hi/lo base

Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides without hip guides	8½–11½ 14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lb)	150