

# 2024 order form for medium

Use dimension chart to select appropriate size.



To order or get a quote:

Go online: <https://www.rifton.com/activity-sizing>

Email: [sales@rifton.com](mailto:sales@rifton.com)

Call: 800.571.8198



		Standard R840 medium		Hi/lo R850 medium							
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$670	<input type="checkbox"/>	R851	\$1390	<input type="checkbox"/>			
		Adjustable w/ back spring	R842	\$760	<input type="checkbox"/>	R852	\$1510	<input type="checkbox"/>			
* Base	CHOOSE	Without spring	R843	\$385	<input type="checkbox"/>	R853	\$2790	<input type="checkbox"/>			
		With spring	R844	\$465	<input type="checkbox"/>	N/A					
* Arm supports (pair)	CHOOSE	Armrests	R868	\$270	<input type="checkbox"/>	R868	\$270	<input type="checkbox"/>			
		Forearm prompts	R869	\$675	<input type="checkbox"/>	R869	\$675	<input type="checkbox"/>			
* Adjustable legs	CHOOSE	Short	R893	\$65	<input type="checkbox"/>	N/A					
		Long	R894	\$85	<input type="checkbox"/>	N/A					
		Short with casters	R895	\$330	<input type="checkbox"/>	N/A					
		Long with casters	R896	\$360	<input type="checkbox"/>	N/A					
* Pads	CHOOSE	Tan	R835	\$275	<input type="checkbox"/>	R835	\$275	<input type="checkbox"/>			
		Blue	R836	\$275	<input type="checkbox"/>	R836	\$275	<input type="checkbox"/>			
		Red	R837	\$275	<input type="checkbox"/>	R837	\$275	<input type="checkbox"/>			
		Pink	R974	\$275	<input type="checkbox"/>	R974	\$275	<input type="checkbox"/>			
		Purple	R975	\$275	<input type="checkbox"/>	R975	\$275	<input type="checkbox"/>			
		Green	R976	\$275	<input type="checkbox"/>	R976	\$275	<input type="checkbox"/>			
Headrest	CHOOSE	Contoured	R897	\$260	<input type="checkbox"/>	R897	\$260	<input type="checkbox"/>			
		Flat	R898	\$260	<input type="checkbox"/>	R898	\$260	<input type="checkbox"/>			
		Adjustable winged	R881	\$290	<input type="checkbox"/>	R881	\$290	<input type="checkbox"/>			
Laterals (pair)	CHOOSE	Small	R828	\$220	<input type="checkbox"/>	R828	\$220	<input type="checkbox"/>			
		Small with chest strap	R829	\$275	<input type="checkbox"/>	R829	\$275	<input type="checkbox"/>			
		Large	R838	\$235	<input type="checkbox"/>	R838	\$235	<input type="checkbox"/>			
Push handles (recommended for all chairs with wheels)			R878	\$166	<input type="checkbox"/>	Included					
Tray (requires armrests)			R668	\$300	<input type="checkbox"/>	R668	\$300	<input type="checkbox"/>			
Handhold			R649	\$31	<input type="checkbox"/>	R649	\$31	<input type="checkbox"/>			
Chest strap (wide)			R863	\$100	<input type="checkbox"/>	R863	\$100	<input type="checkbox"/>			
Butterfly harness	CHOOSE	Regular	R864	\$135	<input type="checkbox"/>	R864	\$135	<input type="checkbox"/>			
		Slim cut	R882	\$135	<input type="checkbox"/>	R882	\$135	<input type="checkbox"/>			
Thigh belt			R865	\$105	<input type="checkbox"/>	R865	\$105	<input type="checkbox"/>			
Pelvic harness			R866	\$190	<input type="checkbox"/>	R866	\$190	<input type="checkbox"/>			
Hip guides (pair)			R867	\$230	<input type="checkbox"/>	R867	\$230	<input type="checkbox"/>			
Abductor	CHOOSE	Small	R884	\$140	<input type="checkbox"/>	R884	\$140	<input type="checkbox"/>			
		Large	R891	\$160	<input type="checkbox"/>	R891	\$160	<input type="checkbox"/>			
Adductors (pair)			R892	\$220	<input type="checkbox"/>	R892	\$220	<input type="checkbox"/>			
Leg prompt			R858	\$245	<input type="checkbox"/>	R858	\$245	<input type="checkbox"/>			
Footboard (recommended for chairs with long legs)			R848	\$525	<input type="checkbox"/>	Included					
Ankle straps (pair)			R899	\$100	<input type="checkbox"/>	R899	\$100	<input type="checkbox"/>			
Sandals (pair)			R954	\$105	<input type="checkbox"/>	R954	\$105	<input type="checkbox"/>			
Wedges (pair, require sandals)			R698	\$80	<input type="checkbox"/>	R698	\$80	<input type="checkbox"/>			
Lumbar and seat support kit			R887	\$88	<input type="checkbox"/>	R887	\$88	<input type="checkbox"/>			
Backrest filler pad			R849	\$135	<input type="checkbox"/>	R849	\$135	<input type="checkbox"/>			
Whitmyer headrest adapter			R804	\$63	<input type="checkbox"/>	R804	\$63	<input type="checkbox"/>			
Directional Assist Wheel retrofit kit			N/A			R965	\$250	<input type="checkbox"/>	R966 (5 kits)	\$1400	<input type="checkbox"/>

\* Required



# Room to grow

## Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"



## Medium chair at its largest dimensions


- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




**Medium standard base**  
R840 Rifton Activity Chair



**Medium hi/lo base**  
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R840 medium Standard base
Frame width	short legs: 23¾
	long legs: 24½
Seat height above floor	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above footboard	short legs: 12½–16½
	long legs: 19–23
Seat angle (tilt-in-space)–degrees	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Backrest angle–degrees	15° forward, 15° back
Footboard knee angle–degrees	5° forward, 20° back
Seat height above footboard	45°–110°
Seat width with hip guides without hip guides	7¾–15¼
Seat depth	8½–11½
Armrest height above seat	14
Distance between laterals	11–16
Backrest height	7–10½
Headrest height above seat	6½–12
Max. working load (lb)	15½–19
	17½–24½
	150

User dimensions (inches)	
Height	40–56
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.</p>	
Item dimensions (inches)	
	R850 medium Hi/lo base
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides without hip guides	8½–11½
Seat depth	14
Armrest height above seat	11–16
Distance between laterals	7–10½
Backrest height	6½–12
Headrest height above seat	15½–19
Max. working load (lb)	17½–24½
	150