

Small Activity Chair 2018



			Standard R820 small		Hi/lo R830 small		
Required components	Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE Adjustable	R821	\$520 <input type="checkbox"/>	R831	\$1,005 <input type="checkbox"/>	
		Adjustable w/ back spring	R822	\$590 <input type="checkbox"/>	R832	\$1,110 <input type="checkbox"/>	
	Base	CHOOSE Without spring	R823	\$305 <input type="checkbox"/>	R833	\$1,810 <input type="checkbox"/>	
		With spring	R824	\$355 <input type="checkbox"/>	N/A		
	Arm supports (pair)	CHOOSE Armrests	R827	\$210 <input type="checkbox"/>	R827	\$210 <input type="checkbox"/>	
		Forearm prompts	R809	\$520 <input type="checkbox"/>	R809	\$520 <input type="checkbox"/>	
	Adjustable legs	CHOOSE	Short	R803	\$41 <input type="checkbox"/>	N/A	
			Long	R894	\$70 <input type="checkbox"/>	N/A	
			Short with casters	R805	\$255 <input type="checkbox"/>	N/A	
			Long with casters	R896	\$285 <input type="checkbox"/>	N/A	
	Pads	CHOOSE	Tan	R815	\$210 <input type="checkbox"/>	R815	\$210 <input type="checkbox"/>
			Blue	R816	\$210 <input type="checkbox"/>	R816	\$210 <input type="checkbox"/>
			Red	R817	\$210 <input type="checkbox"/>	R817	\$210 <input type="checkbox"/>
			Pink	R971	\$210 <input type="checkbox"/>	R971	\$210 <input type="checkbox"/>
Purple			R972	\$210 <input type="checkbox"/>	R972	\$210 <input type="checkbox"/>	
Green			R973	\$210 <input type="checkbox"/>	R973	\$210 <input type="checkbox"/>	
Optional accessories	Headrest	CHOOSE	Contoured	R897	\$210 <input type="checkbox"/>	R897	\$210 <input type="checkbox"/>
			Flat	R898	\$195 <input type="checkbox"/>	R898	\$195 <input type="checkbox"/>
			Adjustable winged	R881	\$260 <input type="checkbox"/>	R881	\$260 <input type="checkbox"/>
	Laterals (pair)	CHOOSE	Small	R828	\$175 <input type="checkbox"/>	R828	\$175 <input type="checkbox"/>
			Small with strap	R819	\$210 <input type="checkbox"/>	R819	\$210 <input type="checkbox"/>
	Push handles (recommended for all chairs with wheels)		R878	\$135 <input type="checkbox"/>	Included		
	Tray (requires armrests)		R667	\$210 <input type="checkbox"/>	R667	\$210 <input type="checkbox"/>	
	Handhold		R649	\$21 <input type="checkbox"/>	R649	\$21 <input type="checkbox"/>	
	Chest strap (wide)		R802	\$64 <input type="checkbox"/>	R802	\$64 <input type="checkbox"/>	
	Butterfly harness		R834	\$100 <input type="checkbox"/>	R834	\$100 <input type="checkbox"/>	
	Thigh belt		R811	\$64 <input type="checkbox"/>	R811	\$64 <input type="checkbox"/>	
	Pelvic harness		R845	\$135 <input type="checkbox"/>	R845	\$135 <input type="checkbox"/>	
	Hip guides (pair)		R814	\$180 <input type="checkbox"/>	R814	\$180 <input type="checkbox"/>	
	Abductor		R884	\$110 <input type="checkbox"/>	R884	\$110 <input type="checkbox"/>	
	Adductors (pair)		R892	\$175 <input type="checkbox"/>	R892	\$175 <input type="checkbox"/>	
	Leg prompt		R854	\$190 <input type="checkbox"/>	R854	\$190 <input type="checkbox"/>	
	Footboard (recommended for chairs with long legs)		R847	\$370 <input type="checkbox"/>	Included		
	Footboard lift		R825	\$200 <input type="checkbox"/>	R825	\$200 <input type="checkbox"/>	
	Mini kit (footboard lift, backrest insert, and trunk support)		R812	\$270 <input type="checkbox"/>	R812	\$270 <input type="checkbox"/>	
	Ankle straps (pair)		R899	\$70 <input type="checkbox"/>	R899	\$70 <input type="checkbox"/>	
	Sandals (pair)		R944	\$77 <input type="checkbox"/>	R944	\$77 <input type="checkbox"/>	
	Wedges (pair, requires sandals)		R697	\$59 <input type="checkbox"/>	R697	\$59 <input type="checkbox"/>	
	Lumbar and seat support kit		R886	\$59 <input type="checkbox"/>	R886	\$59 <input type="checkbox"/>	
	Backrest filler pad		R849	\$110 <input type="checkbox"/>	R849	\$110 <input type="checkbox"/>	
	Whitmyer headrest adapter		R804	\$40 <input type="checkbox"/>	R804	\$40 <input type="checkbox"/>	
	Total			R820	\$ _____	R830	\$ _____

Room to grow

Small chair with mini kit at its smallest dimensions

- 1 Backrest height 12½"
- 2 Trunk support width 5½"
- 3 Seat depth 7"
- 4 Seat width 7"
- 5 Armrest height above seat 5"
- 6 Seat height above footboard 6"




Small chair at its largest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 11½"
- 3 Seat depth 12"
- 4 Seat width 12"
- 5 Armrest height above seat 7½"
- 6 Seat height above footboard 12"




Small standard base
R820 Rifton Activity Chair

User dimensions (inches)	
Height with mini kit	32–48 28–38
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth, and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.	
Item dimensions (inches)	
	R820 small Standard base
Frame width	short legs: 21 long legs: 23
	short legs w/ casters: 22½ long legs w/ casters: 23½
Seat height above floor	short legs: 9½–12½ long legs: 18½–21½
	short legs w/ casters: 13½–16½ long legs w/ casters: 18½–21½
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	9–12 6–9
Seat width with hip guides without hip guides	7–9 12
Seat depth with mini kit	8–12 7–11
Armrest height above seat	5–7½
Trunk support width	5½–11½
Backrest height	12½–15½
Headrest height above seat	14½–21
Max. working load (lbs)	75



Small hi/lo base
R830 Rifton Activity Chair

User dimensions (inches)	
Height with mini kit	32–48 28–38
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth, and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.	
Item dimensions (inches)	
	R830 small Hi/lo base
Frame width	26
Seat height above floor	10–25
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	9–12 6–9
Seat width with hip guides without hip guides	7–9 12
Seat depth with mini kit	8–12 7–11
Armrest height above seat	5–7½
Trunk support width	5½–11½
Backrest height	12½–15½
Headrest height above seat	14½–21
Max. working load (lbs)	75