

# 2024 order form



To order or get a quote:  
 Go online: <https://www.rifton.com/stander-sizing>  
 Email: [sales@rifton.com](mailto:sales@rifton.com)  
 Call: 800.571.8198


Use dimension chart to select appropriate size.

		S410 size 1		S420 size 2		
<b>* Frame</b>		S411	\$1,900	<input type="checkbox"/>	S421 \$2,500 <input type="checkbox"/>	
<b>* Color</b>		<b>CHOOSE</b>	Blue	<input type="checkbox"/>	Blue <input type="checkbox"/>	
			Red	<input type="checkbox"/>	Red <input type="checkbox"/>	
<b>* CHOOSE CONFIGURATION</b>	<b>Prone</b> (with laterals, pelvic harness) HCPCS code: E0638	<b>Prone knee cuffs</b> (pair)		S412 \$500 <input type="checkbox"/>	S422 \$510 <input type="checkbox"/>	
	<b>Supine</b> (with upper trunk support, laterals) HCPCS code: E0638	<b>Head laterals</b> (pair)		S413 \$550 <input type="checkbox"/>	S423 \$560 <input type="checkbox"/>	
		<b>Butterfly harness</b>		S497 \$200 <input type="checkbox"/>	S497 \$200 <input type="checkbox"/>	
		<b>* Hip strap</b>	<b>CHOOSE</b>	Standard	S456 \$300 <input type="checkbox"/>	S457 \$310 <input type="checkbox"/>
				Rotation control	S472 \$300 <input type="checkbox"/>	S473 \$310 <input type="checkbox"/>
		<b>* Knee supports</b>	<b>CHOOSE</b>	Soft straps (pair)	S476 \$300 <input type="checkbox"/>	S477 \$310 <input type="checkbox"/>
	Split supports (pair)			S401 \$200 <input type="checkbox"/>	S402 \$210 <input type="checkbox"/>	
			Split supports (pair)	N/A	S487 \$425 <input type="checkbox"/>	
	<b>Multi-Position</b> (with upper trunk support, laterals, pelvic harness) HCPCS code: E0641	<b>Head laterals</b>		S414 \$850 <input type="checkbox"/>	S424 \$870 <input type="checkbox"/>	
		<b>Butterfly harness</b>		S497 \$200 <input type="checkbox"/>	S497 \$200 <input type="checkbox"/>	
		<b>* Hip strap</b>	<b>CHOOSE</b>	Standard	S456 \$300 <input type="checkbox"/>	S457 \$310 <input type="checkbox"/>
				Rotation control	S472 \$300 <input type="checkbox"/>	S473 \$310 <input type="checkbox"/>
<b>* Knee supports</b>		<b>CHOOSE</b>	Soft straps (pair)	S476 \$300 <input type="checkbox"/>	S477 \$310 <input type="checkbox"/>	
			Split supports (pair)	S401 \$200 <input type="checkbox"/>	S402 \$210 <input type="checkbox"/>	
		Split supports (pair)	N/A	S487 \$425 <input type="checkbox"/>		
<b>* Trunk strap</b>	<b>CHOOSE</b>	Standard	S485 \$200 <input type="checkbox"/>	S485 \$200 <input type="checkbox"/>		
		Rotation control	S462 \$300 <input type="checkbox"/>	S463 \$310 <input type="checkbox"/>		
<b>* Leg type</b>	<b>CHOOSE</b>	Standard	S466 \$300 <input type="checkbox"/>	S467 \$310 <input type="checkbox"/>		
		Advanced	S415 \$650 <input type="checkbox"/>	S425 \$670 <input type="checkbox"/>		
		Advanced	N/A	S426 \$970 <input type="checkbox"/>		
<b>Tray</b>		S491 \$550 <input type="checkbox"/>	S491 \$550 <input type="checkbox"/>			
<b>Additional trunk laterals</b> (pair)		S494 \$200 <input type="checkbox"/>	S494 \$200 <input type="checkbox"/>			
<b>Medial/lateral thigh supports</b> (pair)		S482 \$200 <input type="checkbox"/>	S482 \$200 <input type="checkbox"/>			
<b>Additional medial/lateral thigh supports</b> (pair)		S482 \$200 <input type="checkbox"/>	S482 \$200 <input type="checkbox"/>			
<b>S1 Sandal lift kit</b>		S409 \$50 <input type="checkbox"/>	N/A			
<b>S2 Headrest adapter plate kit</b>		N/A	S427 \$250 <input type="checkbox"/>			

**\* Required**

The HCPCS codes indicated above are suggestions only, based on knowledge of our products and the HCPCS definitions established by CMS. It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. Actual coding and coverage varies by payer; some may accept or require alternative HCPCS codes, including miscellaneous codes, to ensure access for their beneficiaries.

# Stander dimensions

User dimensions (inches)	S410 Size 1	S420 Size 2
Height	28–42*	38–53
 <p><b>Key user dimension: height</b>            Select the appropriate stander by the user's overall height.            Choose the model that allows for growth.            *The smallest users may require the S409 Sandal lift kit for best positioning.</p> <p><b>Important:</b> User's weight must not exceed the maximum working load.</p>		
Item dimensions (inches)	S410 Size 1	S420 Size 2
Base length x width	32½ x 21½	39 x 27
Height of horizontal body support above floor	22–28	22½–36½
Sandal to knee height	6¾–11	10–15¼
Knee to hip height	6¼–10½	8½–13½
Sandal to top of headrest	31½–43	41–55
Width between trunk laterals	5½–10	7–11½
Width between medial/lateral thigh supports	2½–4	3¼–5
Basic item weight PRONE (lb)	47	72
Basic item weight SUPINE (lb)	52	78
Max. working load (lb)	75	100



**Size 1 range of adjustment**



**Size 2 range of adjustment**