

2025 order form

Use dimension chart to select appropriate size.


Scan for sizing and pricing:
<https://www.rifton.com/standers>



		S410 size 1	S420 size 2	S430 size 3	
* Frame		S411 \$1,940 <input type="checkbox"/>	S421 \$2,550 <input type="checkbox"/>	S431 \$3,160 <input type="checkbox"/>	
* Color	CHOOSE	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	
		Red <input type="checkbox"/>	Red <input type="checkbox"/>	Red <input type="checkbox"/>	
* CHOOSE CONFIGURATION	Prone (with pelvic harness)	S412 \$305 <input type="checkbox"/>	S422 \$315 <input type="checkbox"/>	S432 \$325 <input type="checkbox"/>	
	New! Prone upper trunk support (recommended)	N/A	N/A	S408 \$250 <input type="checkbox"/>	
	* Trunk strap	CHOOSE	Standard	S462 \$305 <input type="checkbox"/>	S463 \$315 <input type="checkbox"/>
			Rotation control	S466 \$305 <input type="checkbox"/>	S467 \$315 <input type="checkbox"/>
	Prone knee cuffs (pair)	S485 \$205 <input type="checkbox"/>	S485 \$205 <input type="checkbox"/>	S486 \$225 <input type="checkbox"/>	
	Supine (with upper trunk support)	S413 \$355 <input type="checkbox"/>	S423 \$365 <input type="checkbox"/>	S433 \$375 <input type="checkbox"/>	
	Head laterals (pair)	CHOOSE	Small	S497 \$205 <input type="checkbox"/>	S497 \$205 <input type="checkbox"/>
			Large	N/A	S498 \$215 <input type="checkbox"/>
	* Trunk straps	CHOOSE	Butterfly harness (standard)	S456 \$305 <input type="checkbox"/>	S457 \$315 <input type="checkbox"/>
			Butterfly harness (slim-cut)	N/A	N/A
Trunk strap (standard)			S462 \$305 <input type="checkbox"/>	S463 \$315 <input type="checkbox"/>	
Trunk strap (rotation control)			S466 \$305 <input type="checkbox"/>	S467 \$315 <input type="checkbox"/>	
* Hip strap	CHOOSE	Standard	S472 \$305 <input type="checkbox"/>	S473 \$315 <input type="checkbox"/>	
		Rotation control	S476 \$305 <input type="checkbox"/>	S477 \$315 <input type="checkbox"/>	
* Knee supports	CHOOSE	Soft straps (pair)	S401 \$205 <input type="checkbox"/>	S402 \$215 <input type="checkbox"/>	
		Split supports (pair)	N/A	S487 \$435 <input type="checkbox"/>	
Multi-Position (with head support, pelvic harness)	S414 \$665 <input type="checkbox"/>	S424 \$685 <input type="checkbox"/>	S434 \$705 <input type="checkbox"/>		
Head laterals	CHOOSE	Small	S497 \$205 <input type="checkbox"/>	S497 \$205 <input type="checkbox"/>	
		Large	N/A	S498 \$215 <input type="checkbox"/>	
Butterfly harness	CHOOSE	Standard	S456 \$305 <input type="checkbox"/>	S457 \$315 <input type="checkbox"/>	
		Slim-cut	N/A	N/A	
* Trunk strap	CHOOSE	Standard	S462 \$305 <input type="checkbox"/>	S463 \$315 <input type="checkbox"/>	
		Rotation control	S466 \$305 <input type="checkbox"/>	S467 \$315 <input type="checkbox"/>	
* Hip strap	CHOOSE	Standard	S472 \$305 <input type="checkbox"/>	S473 \$315 <input type="checkbox"/>	
		Rotation control	S476 \$305 <input type="checkbox"/>	S477 \$315 <input type="checkbox"/>	
* Knee supports	CHOOSE	Soft straps (pair)	S401 \$205 <input type="checkbox"/>	S402 \$215 <input type="checkbox"/>	
		Split supports (pair)	N/A	S487 \$435 <input type="checkbox"/>	
New! Prone upper trunk support (recommended)	N/A	N/A	S408 \$250 <input type="checkbox"/>		
Prone knee cuffs (pair)	S485 \$205 <input type="checkbox"/>	S485 \$205 <input type="checkbox"/>	S486 \$225 <input type="checkbox"/>		
* Trunk laterals	CHOOSE	Small	S494 \$205 <input type="checkbox"/>	S494 \$205 <input type="checkbox"/>	
		Large	N/A	S495 \$215 <input type="checkbox"/>	
* Leg type	CHOOSE	Standard	S415 \$665 <input type="checkbox"/>	S425 \$685 <input type="checkbox"/>	
		Advanced	N/A	S426 \$990 <input type="checkbox"/>	
Tray	S491 \$560 <input type="checkbox"/>	S491 \$560 <input type="checkbox"/>	S492 \$610 <input type="checkbox"/>		
Additional trunk laterals	CHOOSE	Small	S494 \$205 <input type="checkbox"/>	S494 \$205 <input type="checkbox"/>	
		Large	N/A	S495 \$215 <input type="checkbox"/>	
Medial/lateral thigh supports (pair)	CHOOSE	Small	S482 \$205 <input type="checkbox"/>	S482 \$205 <input type="checkbox"/>	
		Large	N/A	S483 \$215 <input type="checkbox"/>	
Additional medial/lateral thigh supports (pair)	CHOOSE	Small	S482 \$205 <input type="checkbox"/>	S482 \$205 <input type="checkbox"/>	
		Large	N/A	S483 \$215 <input type="checkbox"/>	
Size 1 sandal lift kit	S409 \$55 <input type="checkbox"/>	NA	NA		
Headrest adapter plate kit	N/A	S427 \$255 <input type="checkbox"/>	S427 \$255 <input type="checkbox"/>		

* Required

Stander dimensions

User dimensions (inches)	S410 Size 1	S420 Size 2	S430 Size 3
Height	28–42*	38–53	48–65
 <p>Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth. *The smallest users may require the S409 Sandal lift kit for best positioning.</p> <p>Important: User's weight must not exceed the maximum working load.</p>			
Item dimensions (inches)	S410 Size 1	S420 Size 2	S430 Size 3
Base length x width	32½ x 21½	39 x 27	43 x 29
Height of horizontal body support above floor	22–28	22½–36½	23½–38
Sandal to knee height	6¾–11	10–15¼	12¼–19
Knee to hip height	6¼–10½	8½–13½	12–17¼
Sandal to top of headrest	31½–43	41–55	49–66½
Width between trunk laterals	5½–10	7–11½	9–14
Width between medial/lateral thigh supports	2½–4	3¼–5	4–5½
Basic item weight PRONE (lb)	47	72	104
Basic item weight SUPINE (lb)	52	78	115
Max. working load (lb)	75	100	185



New! Size 3 Stander at its smallest and largest settings