


Support Station dimensions

User dimensions (inches)		K710
Height		46 – 80
 <p>Important: User height and weight must not exceed the maximum number.</p>		
Item dimensions (inches)		K710
Trunk board size		24 x 24
Maximum trunk board height		45
Minimum trunk board height		26
Height of frame		65½
Width of frame		27½
Depth when trunk board horizontal		31
Depth when collapsed against wall		8
Max. working load (lb)		250

2025 order form

Scan for sizing and pricing:
<https://www.rifton.com/supportstation>



		K710			
* Frame		K715	\$3,890	<input type="checkbox"/>	
Support strap		K716	\$130	<input type="checkbox"/>	
Kneeboard		K717	\$775	<input type="checkbox"/>	
Handholds	CHOOSE	Pair	K718	\$100	<input type="checkbox"/>
		Single	K719	\$50	<input type="checkbox"/>

* Required