

2025 order form

Use dimension chart to select appropriate size.


Scan for sizing and pricing:
<https://www.rifton.com/tram>




		K310 Standard (better for gait training)	K320 Low-base (fits under low beds)	
* Frame includes one battery, one battery charger, accessories tote and charger power cord. (International power cords available on request.)	CHOOSE	Frame with scale	K301 \$7,150 <input type="checkbox"/>	K321 \$7,150 <input type="checkbox"/>
		Frame	K302 \$5,850 <input type="checkbox"/>	K322 \$5,850 <input type="checkbox"/>
* Thigh straps (pair)	CHOOSE	Narrow	K308 \$190 <input type="checkbox"/>	K308 \$190 <input type="checkbox"/>
		Wide	K338 \$190 <input type="checkbox"/>	K338 \$190 <input type="checkbox"/>
Switch pole		K666 \$275 <input type="checkbox"/>	K666 \$275 <input type="checkbox"/>	
Arm supports (pair)	CHOOSE	Small arm prompts	K305 \$580 <input type="checkbox"/>	K305 \$580 <input type="checkbox"/>
		Large arm prompts	K315 \$660 <input type="checkbox"/>	K315 \$660 <input type="checkbox"/>
		Arm platforms	K397 \$435 <input type="checkbox"/>	K397 \$435 <input type="checkbox"/>
		Arm platforms w/ handgrips	K396 \$490 <input type="checkbox"/>	K396 \$490 <input type="checkbox"/>
Pelvic support	CHOOSE	Small	K306 \$295 <input type="checkbox"/>	K306 \$295 <input type="checkbox"/>
		Medium	K316 \$330 <input type="checkbox"/>	K316 \$330 <input type="checkbox"/>
		Large	K326 \$350 <input type="checkbox"/>	K326 \$350 <input type="checkbox"/>
Hip positioner	CHOOSE	Small	K667 \$360 <input type="checkbox"/>	K667 \$360 <input type="checkbox"/>
		Large	K668 \$400 <input type="checkbox"/>	K668 \$400 <input type="checkbox"/>
Hip positioner pad	CHOOSE	Small	K541 \$120 <input type="checkbox"/>	K541 \$120 <input type="checkbox"/>
		Large	K542 \$160 <input type="checkbox"/>	K542 \$160 <input type="checkbox"/>
Additional battery (recommended for institutional use)		K307 \$175 <input type="checkbox"/>	K307 \$175 <input type="checkbox"/>	
Additional thigh straps	CHOOSE	Narrow single	K318 \$95 <input type="checkbox"/>	K318 \$95 <input type="checkbox"/>
		Narrow pair	K308 \$190 <input type="checkbox"/>	K308 \$190 <input type="checkbox"/>
		Wide single	K328 \$95 <input type="checkbox"/>	K328 \$95 <input type="checkbox"/>
		Wide pair	K338 \$190 <input type="checkbox"/>	K338 \$190 <input type="checkbox"/>
Swivel locks (pair)		K319 \$220 <input type="checkbox"/>	K324 \$220 <input type="checkbox"/>	

* Required


TRAM dimensions

User dimensions (inches)		K310 and K320	
Girth		22 min–55 max	
 <p>Key user dimension: girth</p> <p>Important: User's weight or height must not exceed 350 lb or 76 inches</p>			
Item dimensions (inches)	K310 TRAM	K320 TRAM	
Overall length	45½	Same	
Overall width	27½–40	Same	
Overall height	43½–58	Same	
TRAM weight (lb)	70	Same	
Turning diameter	50	Same	
Floor to top of base leg for furniture clearance	6¾	4¼	
Min. user armpit height (when standing or ambulating)	30½	Same	
Max. user height (when standing or ambulating)	76	Same	
Max. working load (lb)	350	Same	

Pelvic support

User dimensions (inches)			
Indicated by the red arrow below			
Key user dimension: half of thigh circumference measured at the groin.			
	K306 small	K316 med.	K326 large
	6–10	8–14	10–18
			
Product dimensions (inches)			
	K306 small	K316 med.	K326 large
a	8½	12	13½
b	16	20	25
c	9	12	14

Hip positioner

Product dimensions (inches)–without pad		
		
	K667 small	K668 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11