

Additional Information for the Rifton Pacer



Chest Prompts

Range of circumference
 small 14"–28"
 med 22"–40"
 large 28"–50"

Height
 small 4½"
 med 6½"
 large 8"



Arm Prompts

Padded surface area
 small 4½" x 7"
 large 5" x 8½"

Adjustment range
 (elbow to fist): small 8"–12½"
 large 10½"–16"



Arm Platforms

Padded surface dimensions: 13½" x 5½"



Thigh Prompts

Length (from clamp to strap)
 small 5"–8"
 large 5"–11"

Max. leg circumference
 small 15"
 large 20"



Communication tray

Inside dimensions: 8½" x 11"
Max. working load: 10 lb.


Pelvic Support



User dimensions (inches)
 Indicated by the red arrow below

Key user dimension: half of thigh circumference.

	K556 small	K557 med.	K558 large
	4–8	6–10	8–14




Product dimensions (inches)

	K556 small	K557 med.	K558 large
a	6½	8½	12
b	12	16	20
c	7	9	12

Hip Positioner



Dimensions (inches)–without pad



	K513 small	K533 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11