



Adaptive Desk

The Rifton Adaptive Desk is a height-adjustable, wheelchair-accessible desk for inclusive classrooms and therapy spaces.



Scan here to order, request a quote, or see detailed product information.

- Supports students of all abilities with comfort, safety, and ease
- Promotes participation, posture, and independence across learning environments
- Thoughtfully designed for everyday use



Key features

1 Inclusive and ADA-compliant

The 40" x 24" desk (R250) features a wide, open base and meets ADA requirements—ideal for most wheelchairs and Rifton seating. The 32" x 24" desk (R240) offers the same function in a narrower size. Check sizing to ensure fit.



2 Therapeutic positioning

The sit-to-stand height range supports seated and standing work to encourage better posture and engagement.



The Adaptive Desk in use



The wide frame accommodates power wheelchairs and large Rifton Activity Chairs.



Tool-free height adjustments make quick sit-to-stand transitions easy.



Scan for sizing
and pricing.

3 Sized to fit

Choose a narrow top for tight spaces, or a wide top to accommodate larger wheelchair and seating systems.



4 Surface tilt

The optional 0–50° surface tilt adjusts with a simple hand lever.



The compact 32" x 24" top fits efficiently in limited classroom space.



The surface tilt supports comfortable access to materials and devices in seated or standing positions.

Adaptive Desk options

Surface width options

Choose between a narrow 32" x 24" surface width and a wider 40" x 24" surface width.



32" x 24" Adaptive Desk R240

The compact top fits efficiently in limited classroom space.

Height range: 23–43"



40" x 24" Adaptive Desk R250

The wider top accommodates most wheelchairs and larger seating systems.

Height range: 23–43"

Base options



Mobile base


The locking casters provide mobility when you need it and stability when you don't.



Stationary base

The stationary feet create a stable platform and match the height of the mobile base.

Surface options



Maple

Black

Gray

Durable, cleanable work surfaces, these raised edge tops are available in maple, black, gray, or magnet-ready black surface. Plain edge tops are available in maple or gray.

Tilt option



The optional 0–50° continuous surface tilt supports reading, writing, and device use in seated and standing positions.

Adaptive Desk components



Device mount adapter

Movable attachment point for communication devices on raised edge desks.



Slide-on ledge

Pairs with the tilt feature to help keep items from sliding when the surface is angled.



Anti-tip rollers

Prevent the desk from tipping forward.



Side storage box

Convenient, within-reach accessory storage.



Elbow supports

Adjustable forearm supports for customized positioning, more versatile than fixed desk cutouts.



Flexible dry erase magnetic sheet

Unroll on the magnet-ready surface for drawing with dry-erase markers; rolls up for storage.



Multifunction magnetic ledge

Attaches anywhere on the magnet-ready surface to hold a device, papers, or visuals in place.

Need support?

The Rifton website offers everything you need:

- Evidence-based articles
- Sample letters of medical necessity
- Video tutorials and training resources
- Product instruction guides
- Clinical resources and best practices



Scan here to order, request a quote, or see detailed product information.