



Mobile Stander

Experience mobility, independence, and exploration with the **Rifton Mobile Stander**. While users discover their world, they're building essential motor skills that support progression toward independent standing.



Scan here to order,
request a quote, or
see detailed product
information.



Key features

1 Progressive weight bearing

The adjustable seat pad enables gradual weight-bearing increases while building essential standing mobility skills.



2 Trunk development

The adjustable body support reduces assistance when lowered as users develop independent trunk control.



The Mobile Stander in use



Large wheels allow users to independently explore their environment.



Position users in prone or supine.



Scan for sizing
and pricing.

3 Social interaction

Enable natural peer interaction and environmental exploration while developing weight-shifting, balance, and postural control.



4 Independent mobility

Optional large wheels attach easily, giving users control and freedom to explore independently.



The additional seat pad enhances back support.



Experience the freedom of standing.

Mobile Stander options



Scan for sizing
and pricing.

Frame size options



Large K170



Medium K150



Small K130



Mini K110

Frame color options



Blue



Lime



Pink



Mobile Stander components



Large wheels

Wheels attach and remove with one-button operation. Wheel locks ensure security, spoke guards protect fingers, and hand rims keep hands clean.



Mounting bracket

Secures arm prompts, hand loops, or a communication tray.



Arm prompts

Provide fully adjustable support for users needing head and shoulder stabilization. (Require mounting bracket.)



Hand loops

Provide secure, adjustable support that enhances user stability and confidence during movement. (Require mounting bracket.)



Communication tray

Adjusts for perfect device positioning and easy user access. (Requires mounting bracket.)



Additional seat pad

Provides enhanced back support for improved comfort. Available in two sizes.



Sandals and wedges

Adjustable sandals position feet securely. Optional wedges provide heel or toe elevation up to 1¾".



Conversion kit

Conversion kit upgrades the K110 mini stander to K130 small size for less than 25% of a new stander's cost.

Need support?

The Rifton website offers video tutorials, training resources, sample letters of medical necessity, evidence-based articles, and more.



Scan here to order, request a quote, or see detailed product information.