

# Rifton Activity Chair

A perfect complement to your feeding therapy practice





**The Rifton Activity Chair**, the first choice of feeding therapists across the country. It's easily configurable with a wide range of support options to help you position your clients with all their unique positioning challenges.



Hi/lo base Activity Chair (R850)



Standard base Activity Chair (R820)

**1 Headrest adjustability** can make or break your feeding success - we offer three styles and even a Whitmyer adapter for compatibility with most Whitmyer headrest options.

**2 Wide range of components** including various pelvic and trunk supports to meet the requirements of every individual.

**3 The large tray** snaps into place with one hand. With a cutout for a snug fit around the user, the tray offers arm support as well as room for mealtime items.

**4 Seat and backrest padding** covered with polyurethane-coated fabric in six color choices is easy to clean and removable for scrubbing.

**5 Hi/lo adjustment** allows for positioning low to the floor and up at the table with family or peers, while the user is in the chair. Tilt-in-space quickly adjusts to the desired angle.

► **Only pay for the components your client needs. Rifton's customer service team is ready to answer questions regarding the available options.**



### More features of the Rifton Activity Chair:

- **Easy forward positioning** encourages active postural participation.
- Ease of **posterior tilt-in-space** or recline for relief of fatigue, for medical respiratory issues, and for a restful digestion period.
- **Tool-free adjustments** make positioning easy, even while the user is in the chair.
- **Range of seat-to-back angles** facilitate trunk and head control and reduce tone.

► Therapists, contact us for a quote and to see if you qualify for an institutional discount.

800.571.8198 | [www.rifton.com](http://www.rifton.com) | [sales@rifton.com](mailto:sales@rifton.com)



"I love that this chair has so much growth potential; it will last a long time for most children. I love that you can measure the angle of the tilt-in-space, and the backrest recline. This improves accuracy of carry-over to home caregivers. The chair has been highly recommended for children with feeding difficulties. It provides the support that is missing in standard high chairs, booster and toddler seats. Appropriate positioning is a key element for success with feeding."

**Carol Elliott, OTR/L**  
**Registered Occupational Therapist**



**Tip:** A well-written Letter of Medical Necessity will assist in obtaining coverage. Download a sample here:

<https://www.rifton.com/feedingchairlmm>