



Support Station

Developed with therapists, **The Rifton Support Station** enables clients to take part in their own hygiene, promoting dignity and independence.



Scan here to order, request a quote, or see detailed product information.

- Increases caregiver safety
- Builds client strength and skills
- Quicker than horizontal diapering



The Support Station folds flat against the wall when not in use.



Key features


1 Dignity and independence

The Support Station empowers individuals with limited weight-bearing to actively participate in their own care, restoring dignity and confidence.



2 Improved caregiver safety

By reducing heavy lifting and awkward positioning, caregivers stay safer, while clothing and skin are easier to access.



Video: See how the Support Station works.

The Support Station in use



Vertical positioning for hygiene promotes dignity and independence by enabling clients to participate in their own care.



The trunk board adjusts to promote a more upright posture and increased weight bearing.



Scan for sizing
and pricing.

3 Strength and skill building

The Support Station encourages clients to use their own strength in transfers, building weight-bearing ability and motor skills over time.



4 Pivot or fixed configurations

Use the pivot configuration (no kneeboard) for rotating transfers to a fixed toilet or commode. Use the fixed configuration (with kneeboard) for upright hygiene tasks and transfers to rolling equipment.



Handholds and slots on the trunk board let users assist with transfers.



The pivot configuration (shown) enables rotating transfers from a wheelchair to a fixed toilet for users with adequate lower-body strength.