

# Product Comparison: Large Pacer, XL Pacer, E-Pacer



**Large Pacer**

**XL Pacer**

**E-Pacer**

|  |         |         |         |
|--|---------|---------|---------|
| <b>Weight limit</b>  | 200 lbs | 250 lbs | 350 lbs |
| <b>Width of base</b><br>(to determine min. door width)                     | 28"     | 31 ½"   | 31 ½"   |
| <b>Adjustment range</b><br>(floor to top of arm prompt)                    | 32"–47" | 34"–49" | 24"–61" |
| <b>Zero-lift sit-to stand and seated transfers</b> (powered height adjust) |         |         | ✓       |
| <b>Bluetooth-equipped scale</b>  |         |         | ✓       |
| <b>Expandable base frame</b>   |         |         | ✓       |
| <b>Caregiver Push Handle (optional)</b>                                    |         |         | ✓       |
| <b>Collapsible frame</b>   | ✓       | ✓       |         |
| <b>Dynamic weight shifting</b>   | ✓       | ✓       |         |
| <b>Thigh prompts</b>   | ✓       | ✓       |         |
| <b>Multi-position saddle (MPS)</b>   | ✓       | ✓       |         |
| <b>Odometer</b>  | ✓       | ✓       | ✓       |
| <b>Ankle prompts</b>   | ✓       | ✓       | ✓       |
| <b>Four-function casters</b>   | ✓       | ✓       | ✓       |
| <b>Hip positioner</b>  | ✓       | ✓       | ✓       |
| <b>Pelvic support</b>  | ✓       | ✓       | ✓       |
| <b>Arm prompts</b>   | ✓       | ✓       | ✓       |